

# Chem competent cells

---

## Introduction

Get started by giving your protocol a name and editing this introduction.

## Materials

- › *E. coli* strain
- › LB medium
- › 0.1 M CaCl<sub>2</sub> solution (ice cold)
- › LB plates with proper antibiotic
- › 0.1 M CaCl<sub>2</sub> solution containing 15% glycerol (ice cold)
- ›

## Procedure

### The day before

1. Note: Make 1/2 portion (amounts bolded)
2. Put the 0.1 M CaCl<sub>2</sub> solution and 0.1 M CaCl<sub>2</sub> solution containing 15% glycerol at 4 °C.
3. Inoculate one single colony of *E. coli* strain in 5 (**2,5**) mL LB medium. Shake at 37 °C overnight.
4. Put at least thirty 1.5 mL Eppendorf tubes at -80 °C.
5. Inoculate 1 (**0,5**) mL overnight culture in 100 (**50**) mL LB medium within a 500 mL flask.
6. Subculture at 37 °C with shaking till OD<sub>600</sub> reaches ~ 0.25-0.3 (about 2 hours subculture time). First check 1,5 h after the inoculation.
7. Chill the culture on ice for 15 minutes.
8. Separate 100 mL chilled bacterial culture into two 50 mL Falcon tubes and centrifuge at 4 °C at 4000 rpm for 10 minutes.
9. Discard the supernatant and resuspend the pellet with 40 (**20**) mL ice-cold 0.1 M CaCl<sub>2</sub> solution.
10. Keep cells on ice again for 30 minutes.
11. Centrifuge cells at 4000 rpm at 4 °C for 10 minutes.
12. Discard the supernatant and resuspend pellet with 5 (**2,5**) mL ice-cold 0.1 M CaCl<sub>2</sub> solution containing 15% glycerol.
13. Pipet 50 µL of cell suspension into -80 °C frozen eppendorfs and directly transfer them to -80 °C freezer.
14. <http://www.genomearchitecture.com/protocols/Ecoli-heat-shock-transformation.html>