

**iGEM TU/e 2015**

Biomedical Engineering

Eindhoven University of Technology  
Room: Ceres 0.04  
Den Dolech 2, 5612 AZ Eindhoven  
The Netherlands  
Tel. no. +31 50 247 55 59  
[2015.igem.org/Team:TU\\_Eindhoven](http://2015.igem.org/Team:TU_Eindhoven)

**Date**

11 August 2015

# **Interlab study Protein Expression Measurement**

## Table of contents

Title	1	Protein Expression Measurement	3
Notebook	1.1	Measurement of the OD600	Error! Bookm
	1.2	Measurements with the plate-reader	Error! Bookm


# 1 Protein Expression

## 1.1 Measurement of the OD600

- 1) Fill a 1.5 mL cuvet with 1mL LB and perform a blanc measurement at OD600.
- 2) Insert 1 mL bacterial culture into another cuvet.
- 3) Measure the OD at 600 nm.
- 4) Calculate the correct dilution to fill the well plates with an OD of 0.5.
  - $X = \text{Measured OD} / 0.5$
  - $1 / X = \text{Volume of added bacterial culture (mL)}$
  - $\text{LB (pH 7.2) added (mL)} = 1 - \text{Volume of added bacterial culture}$
- 5) Fill a new 1.5 mL cuvet with the bacterial culture and LB (pH 7.2) in the calculated proportions.
- 6) Measure the OD at 600 nm again.
- 7) Add up LB (pH 7.2) until an OD of 0.5 ( $\pm 0.025$ ) is reached.

## 1.2 Measurements with the plate-reader

- 1) Start the Tecan infinite F500 and heat up to 37°C.
- 2) Fill row A of a 96 black well plate flat bottom with 200µL LB (pH 7.2).
- 3) Fill the plate with 200 µL of the diluted constructs and the various controls.
  - In our measurements the lay-out was as follows:
    - Row A: LB medium
    - Row B: Colony 1,2 and 3 of construct 1
    - Row C: Colony 1,2 and 3 of construct 2
    - Row D: Colony 1,2 and 3 of construct 3
    - Row E: Colony 1, 2 and 3 of construct 4: Negative control – Empty BL21 bacteria
    - Row F: Colony 1, 2 and 3 of construct 5: Positive control – GFP with working promotor.
    - Row G: Colony 1, 2 and 3 of construct 6: Negative control – BL21 without GFP and with B0062 promotor.



- 4) Measure the plate with the following settings:
  - a. Excitation filter: 485nm (20)
  - b. Emission filter: 535nm (25)
  - c. Heated to 37°C
  - d. Orbital shaking at 2.5mm amplitude and frequency 247.6 orbital for 180 seconds.
  - e. GFP gain 30
- 5) Safe the results