

Gibson assembly protocol

1. Measure the DNA concentration (ng/ μl) of each assembly piece.
2. Put the 15 μl Gibson mix on ice.
3. Transfer units to pmol/μl by the following formula:

$$DNA\ Con' \left[\frac{Pmol}{\mu l} \right] = \frac{DNA\ Con' \left[\frac{ng}{\mu l} \right] \times 1000}{Part\ size \times 650}$$

4. Set up the following reaction on ice:

התוצאות של המבחן הן כאלו שציינתם. המבחן הוא על ידי חישוב המוליות של ה-DNA. המוליות היא כמות המולקולות של ה-DNA. המוליות היא כמות המולקולות של ה-DNA. המוליות היא כמות המולקולות של ה-DNA.

	Recommended Amount of Fragments Used for Assembly		
	2-3 Fragment Assembly	4-6 Fragment Assembly	Positive Control**
Total Amount of Fragments	0.02-0.2 pmols* X μl	0.2-1 pmols* X μl	10 μl
Gibson Assembly Master Mix (2X)	15 10 μl	15 10 μl	15 10 μl
Deionized H ₂ O	10-X μl	10-X μl	0
Total Volume	20 μl***	20 μl***	20 μl

* 50 ng of 5,000 bp dsDNA is about 0.015 pmols.
 50 ng of 500 bp is about 0.15 pmols.
 Optimized cloning efficiency is 50-100 ng of vectors with 2-3 fold of excess inserts.
 Use 5 times more of inserts if size is less than 200 bps.
 ** Control reagents are provided for two experiments.
 *** If greater numbers of fragments are assembled, additional Gibson Assembly Master Mix may be required.

*Use 0.05 pmol/μl of the plasmid and 0.15 pmol/μl of the insert

5. Incubate samples in a thermocycler at 50 °C for 60 minutes.
 Following incubation, store samples on ice or at -20 °C for subsequent transformation.
6. Remove 2-5 μl of the assembly product and transform into competent cells of interest.