

Gene synthesis using PCR amplification of primers

Reagents

- 28uL ddH₂O
- 10uL S₁₅ buffer (5x): HuTao
- 1uL dNTP mix (10mM mix) (Synbio Tech Co., Ltd.)
- Primers (20pM stock)
- 1uL S₁₅TM DNA polymerase (2.5U/μl) (Synbio Tech Co., Ltd.)

A. PCR I.

1. Take 3uL of each primer into a 1.5mL centrifuge tube
2. Add ddH₂O into the tube (to make up the volume to 200uL)
3. Finger flick to mix, then centrifuge briefly.
4. Prepare the PCR system as the following order:
28uL ddH₂O
10uL S₁₅ buffer (5x)
1uL dNTP mix (10mM mix)
10uL primer mix (20pM stock)
1uL S₁₅TM DNA polymerase (2.5U/μl)
5. Finger flick to mix, then centrifuge briefly.
6. Place the reaction tubes in the Thermal Cycler.
7. PCR will be run with the following programme:
1 cycle of 98 °C, 2 min
18cycles of 98 °C, 10 sec
 58 °C, 20 sec
 72 °C, 30 sec
Final extension 72 °C, 1 min
Hold at 4 °C

B. PCR II

1. Prepare the PCR system as the following order:
35uL ddH₂O
10uL S₁₅ buffer (5x)
1uL dNTP mix (10mM mix)
1uL forward primer the first (20pM stock)
1uL reverse primer the last (20pM stock)
1uL S₁₅TM DNA polymerase (2.5U/μl)
2. Finger flick to mix, then centrifuge briefly.
3. Place the reaction tubes in the Thermal Cycler.
4. PCR will be run with the following programme:
1 cycle of 98 °C, 2 min
25cycles of 98 °C, 10 sec
 58 °C, 20 sec
 72 °C, 30 sec
Final extension 72 °C, 1 min

Hold at 4 °C

5. Run the detection gel electrophoresis to test the results.