Protocol for RNA synthesis with NEB HiScribe T7 Quick High

Yield RNA Synthesis Kit

Thaw the necessary kit components, mix and pulse-spin in microfuge to collect solution to

the bottoms of tubes. Keep on ice.

Assemble the reaction at room temperature in the following order:

Nuclease-free water X μI

NTP Buffer Mix 10 µl (6.7 mM each NTP final)

Template DNA $X \mu I (1 \mu g)$

T7 RNA Polymerase Mix 2µl

Total reaction volume 30µl

Mix thoroughly and pulse-spin in a microfuge. Incubate at 37°C for 2 hours.

To remove template DNA, add 30 µl nuclease-free water to each 20µl reaction, followed

by 2µl of DNasel (RNase-free), mix and incubate for 15 minutes at 37°C.

Reference: NEB HiScribe T7 Quick High Yield RNA Synthesis Handbook