

iGEM TU/e 2015

Biomedical Engineering

Eindhoven University of Technology

Room: Ceres 0.04

Den Dolech 2, 5612 AZ Eindhoven

The Netherlands

Tel. no. +31 50 247 55 59

2015.igem.org/Team:TU_Eindhoven

Preparing FACS samples

Table of contents

Preparing FACS samples	1	Preparing FACS samples	3
	1.1	Materials	3
	1.2	Setup & Protocols	3

1 Preparing FACS samples

Estimated bench time: 90 minutes

Estimated total time: 3 hours

Purpose: Preparing the bacteria (after protein expression) for a FACS measurement. The bacteria will be extracted from the culture media and a fluorescent dye is added to covalently bind to the proteins.

1.1 Materials

- 1,5 ml cuvettes
- 5 mM DBCO-PEG₄-TAMRA
- Cell Density Meter (OD600)
- ddH₂O
- Eppendorf tubes
- MiniSpin centrifuge
- PBS-0.5%BSA
- Pipettes and tips
- Shaking block
- Tabletop centrifuge
- Tin foil

1.2 Setup & Protocols

- Spin down the cells in the tabletop centrifuge for 15 minutes at 3,000 xg and 4°C.
- Discard supernatant.
- Resuspend with 1 ml PBS-0.5%BSA and transfer to a 1.5 ml Eppendorf tube.
- Spin down the cells in the MiniSpin centrifuge for 1 minute at 13,400 rpm.
- Discard supernatant
- Resuspend with 1 ml PBS-0.5%BSA.
- Perform an OD measurement on a 20x dilution of the culture sample.
OD measurement (OD600):
 - Blank: 950 µl ddH₂O and 50 µl PBS-0.5%BSA
 - Sample: 950 µl ddH₂O and 50 µl PBS-0.5%BSA
- Multiply the OD with 20.
- Calculate the amount of cells in the culture using the Agilent Technologies website.¹
- Make a dilution with a concentration of 1×10^9 cells/ml.
- Prepare tubes for the FACS (mix well & cover the samples).

Tube	[DBCO]	Cells [1×10^9] (µl)	DBCO-TAMRA (5mM) (µl)	Contains pAzf
1	0	200	-	Yes
2	30 µM	200	1.21	Yes
3	30 µM	200	1.21	No

- Incubate the samples in a shaking block for 1 hour at 300 rpm and 4°C. Make sure the tubes are in the dark.
- Spin down the cells in the MiniSpin centrifuge for 10 minutes at 13,400 rpm.

¹ http://www.genomics.agilent.com/biocalculators/calcODBacterial.jsp?_requestid=826255

- Discard the supernatant.
- Resuspend with 1 ml ice-cold-PBS-0.5%BSA. (keep the tubes on ice)
- Spin down the cells for 10 minutes at 13,400 rpm.
- Discard the supernatant and put the pellet on ice in the dark until FACS.
- For the FACS: resuspend with 200 μ l ice-cold PBS-0.5%BSA.