

3 α -HSD Activity Measurement Protocol

Goal: verification of enzymatic activity of 3 α -HSD by measurement of decrease in NADPH fluorescence over time in presence of dihydrotestosterone (DHT).

Protocol:

Preparation:

1. **The day before:** prepare starters in LB (5ml with 5 μ l antibiotics) and place in shaker overnight at 37°C (16-20hr)
2. Dilute the bacteria 1/100 (50ml LB + 50 μ l antibiotics) in a 250ml Erlenmeyer and check O.D. at 600 nm.
3. Grow culture until it reaches O.D.(600nm) = 0.6.
4. Transfer to 50ml falcon and centrifuge for 5 min at 5000rpm. Re-suspend the pellet in BA (bioassay) medium.
5. Add 0.1mM IPTG (5 μ l) and place in shaker for 2-3 hr.
6. Transfer 1ml into Eppendorf tube and centrifuge for 5 min at 5000rpm.
7. Re-suspend the pellet in 1ml phosphate buffer (PBS) pH=7.4.
8. Sonicate at 20% amplitude for 3 cycles of: 5 sec sonication and 30 sec on ice.
9. Centrifuge for 10 minute at max speed.

Reaction:

1. Keep 96-well microplate on ice and add 170 μ l of the cell lysate supernatant.
2. Add 150 μ M NADPH to each plate as fast as possible to minimize degradation.
3. Place plate in plate reader for 30min at 37°C for stabilization. Use kinetic measurement program to measure fluorescence every 1 minute.
4. Take out plate and add 50 μ M DHT to each plate using multi-pipette.
5. Place plate in plate reader for 3hr at 37°C and measure every 1 minute.

Comments:

1. For NADPH We used 1.5mM stock solution, made by dilution of powder purchased from Sigma-Aldrich in 0.01N NaOH. Link to product:
<https://www.sigmaaldrich.com/catalog/product/sigma/n1630?lang=en®ion=IL>
2. For DHT we used 1mM stock solution, made by dilution of 1mg/ml solution purchased from Sigma-Aldrich in methanol. Link to product:
<http://www.sigmaaldrich.com/catalog/product/cerillian/d073?lang=en®ion=IL>

3. Fluorescence measurement was performed in Tecan-infinite 200Pro plate reader, using 340nm excitation and 460nm emission filter.
4. Controls: no 3αHSD gene, no lysate, no DHT, no NADPH.

Written by Ruth, Alexey and Adi