



# Protocols and Procedures

NAIT iGem 2015

# Table of Contents

<b>PCR</b> .....	<b>2</b>
Make PCR Master Mix .....	2
Set Up Individual Reactions.....	2
Run PCR .....	2
Run Agarose gel electrophoresis to confirm that PCR was successful .	3
QIAquick PCR Purification .....	3
<b>Digestion and Ligation</b> .....	<b>4</b>
Digestion of vectors and inserts .....	4
rSAP Protocol.....	4
Ligation of Vectors and Inserts .....	4
<b>Transforming Bacteria</b> .....	<b>5</b>
Transformation .....	5
<b>Validating the Transformation</b> .....	<b>6</b>
Colony PCR (KAPA HiFi Hotstart) .....	6
Agarose Gel Protocol (1 Gel).....	6
Sample Prep .....	6
Review Results .....	7
<b>Protein Expression, Isolation and Purification</b> .....	<b>7</b>
Protein expression using BL21(DE3) .....	7
Protein Purification under Denaturing .....	7
<b>SDS PAGE and Silver Staining</b> .....	<b>8</b>
Make PA Gel .....	8
Run Electrophoresis.....	9
Silver Stain Plus.....	10

# PCR

## Make PCR Master Mix

Before beginning, ensure all reaction components are properly thawed and mixed.

- I. Calculate the required volumes of each component based on the following table:

Component	50 $\mu$ L Reaction	Final Concentration
PCR Grade Water	Up to 50 $\mu$ L	N/A
2X KAPA HiFi HS	25.0 $\mu$ L	1X
10 $\mu$ M Forward Primer	1.5 $\mu$ L	0.3 $\mu$ L
10 $\mu$ M Reverse Primer	1.5 $\mu$ L	0.3 $\mu$ L
Template DNA	As required	As Required

- Reaction Volumes may be adjusted between 10 - 50  $\mu$ L

- < 1 ng less complex DNA (0.1 to 1.0 ng) per 50  $\mu$ L

## Set Up Individual Reactions

- I. Transfer the appropriate volumes of PCR master mix, template, and primer to individual PCR tubes
- II. Cap or seal individual reactions.
- III. Mix and centrifuge briefly.

## Run PCR

**NOTE:** A PCR gradient was conducted to determine the optimal annealing temperature needed for best results. We observed that the best annealing temperature was 61°C.

- I. Perform PCR with the following Cycle Protocol:

Step	Temperature	Duration	Cycles
Initial Denaturation	95°C	3 min	1
Denaturing	98°C	20 sec	25
Annealing	61°C	30 sec	25
Extension	72°C	15 sec	25
Final Extension	72°C	1 min	1

**NOTE:** PCR products can be left overnight at 4°C

## Run Agarose gel electrophoresis to confirm that PCR was successful

### Agarose gel Protocol

- I. Mix 99mL of TA buffer with 1g of Agarose and mix on a hotplate until boil
- II. Add 10 $\mu$ L of ethidium bromide after boil and pour into Agarose casting tray
- III. Insert comb at desired depth
- IV. Let stand until gel is formed
- V. Add samples into well of desire on Agarose gel
- VI. Run gel at 70V for desired time

### QIAquick PCR Purification

- I. Add 5X Buffer PB to 1X of the PCR reaction and mix into a separate 1.5 mL micro centrifuge tube. If the colour of the mixture is orange or violet, add 10  $\mu$ K 3M sodium acetate, pH 5.0 and mix. The colour of the mixture will turn yellow
- II. Place a QIAquick column into a provided 2 mL collection tube and then place into centrifuge and transfer mix
- III. To bind DNA, apply the samples to the QIAquick column and centrifuge at 13 000 rpm from 60s. Discard the fluid that is in the collection tube (aka the flow-through). Place column back into the same tube.
- IV. Centrifuge the column once more in the provided 2 mL collection tube for 1 minute to remove residual wash buffer
- V. Dab the bottom of the column on Kim Wipe tissue.
- VI. Place each column in a clean 1.5 mL micro centrifuge tube
- VII. Add 50  $\mu$ L of EB to the column (Aim for the center of the membrane inside the column)
- VIII. Centrifuge at 13 000 RPM for 60 s
- IX. DO NOT THROW OUT FLOW THROUGH. Using a micropipette, set at 65  $\mu$ L, take flow through and add it once more to the column
- X. Centrifuge again (13 000 RPM, 60s)
- XI. Dispose of column (Purified DNA should now be in the 1.5 mL microfuge tube)

## Digestion and Ligation

### Digestion of vectors and inserts

Two enzymes are needed in order to cut the DNA stands at specific restriction sites

- I. Find the protocol of restriction enzymes being used
- II. Mix the following components:

Enzyme 1	1 $\mu$ L
Enzyme 2	1 $\mu$ L
DNA	1 $\mu$ g
Buffer	5 $\mu$ L
Nuclease Free H <sub>2</sub> O	Up to 50 $\mu$ L

\*Ratio of enzymes to DNA is 1:1 and water is added to maintain the concentration based on each enzyme

- III. Let reaction occur for 1 hour

### rSAP Protocol

Removal of phosphorylated ends

Preformed after digestion of psB1C3

#### Protocol

- I. Prepare sample using the following

rSAP	1 Unit
DNA	1pmol of DNA ends
Nuclease Free H <sub>2</sub> O	Up to 20 $\mu$ L
CutSmart Buffer (10X)*	2 $\mu$ L

\*Provided with rSAP

[total volume should equal 20 $\mu$ L]

- II. Incubate at 37°C for 30 minutes
- III. Stop the reaction by heat inactivation at 65°C for 5 minutes

### Ligation of Vectors and Inserts

- I. Set up the following reaction mixture on ice:

Component	20 $\mu$ L Reaction
10X T4 Ligase Buffer	2 $\mu$ L
Vector DNA (4kb)	50 ng (0.020 pmol)

Insert DNA (1kb)	37.5 ng (0.60 pmol)
Nuclease Free H <sub>2</sub> O	Up to 20 $\mu$ L
T4 DNA Ligase	1 $\mu$ L

- II. Gently mix the reaction by pipetting up and down and microfuge briefly
- III. For cohesive ends, incubate at 16°C overnight or at room temperature for 10 minutes
- IV. For blunt ends or single base overhangs, incubate at 16°C overnight or room temperature for 2 hours
- V. Heat inactivate at 65°C for 10 minutes
- VI. Chill on ice and transform 1-5  $\mu$ L of the reaction into 50  $\mu$ L competent Cells

## Transforming Bacteria

### Transformation

- I. Thaw tube of BL21 (DE3) competent E. coli cells on ice for 10 min
- II. Add 1-5  $\mu$ L containing 1pg-100ng of plasmid DNA to the cell mixture. Carefully flick the tube 4-5 times to mix cells and DNA. **(DO NOT VORTEX)**
- III. Place the mixture on ice for 30 minutes **(DO NOT MIX)**
- IV. Heat shock at exactly 42°C for exactly 10 seconds. **(DO NOT MIX)**
- V. Place on ice for 5 Minutes. **(DO NOT MIX)**
- VI. Pipette 950 $\mu$ L of room temperature SOC into the mixture.
- VII. Place at 37°C for 60 minutes. Shake vigorously (250RPM) or rotate
- VIII. Warm selection plates to 37°C
- IX. Mix the cells thoroughly by flicking the tube and inverting, then perform several 10-fold serial dilutions in SOC
- X. Spread 50-100 $\mu$ L of each dilution onto a selection plate and incubate overnight at 37°C. alternatively, incubate at 30°C for 24-36 hours or 25°C for 48 hours

# Validating the Transformation

## Colony PCR (KAPA HiFi Hotstart)

- I. Prepare the PCR master mix using the following

Component	50µL Reaction*	Final Concentration
PCR Grade Water	Up to 50 µL	N/A
2X KAPA HiFi HS	25.0 µL	1X
10 µM Forward Primer **	1.5 µL	0.3 µL
10 µM Reverse Primer**	1.5 µL	0.3µL
Colony DNA	As required	As Required

\* Reaction Volumes may be adjusted between 10 - 50 µL

\*\* Provided by iGem

- II. Transfer the appropriate volumes of PCR master mix and colony DNA to individual PCR tubes
- III. Cap or seal individual reactions.
- IV. Mix and centrifuge briefly.
- V. Perform PCR using the following cycle protocol

Step	Temperature	Duration	Cycles
Initial Denaturation	95°C	3 min	1
Denaturing	98°C	20 sec	25
Annealing	61°C	30 sec	25
Extension	72°C	15 sec	25
Final Extension	72°C	1 min	1

## Agarose Gel Protocol (1 Gel)

- I. Mix 99mL of TA buffer with 1g of Agarose and mix on a hotplate until boil
- II. Add 10µL of ethidium bromide after boil and pour into Agarose casting tray
- III. Insert comb at desired depth
- IV. Let stand until gel is formed
- V. Add samples into well of desire on Agarose gel
- VI. Run gel at 70V for desired time

## Sample Prep

- I. Add 10µL sample with 2µL sample buffer for each well. ~1:5 ratio

**NOTE:** allow for excess reagents if needed

## Review Results

- I. Any samples under 339 base pairs in length were determined to be unsuccessful as they had no insert

## Protein Expression, Isolation and Purification

### Protein expression using BL21(DE3)

- I. Transform expression plasmid into BL21 (DE3). Plate on antibiotic selection plates and incubate overnight at 37°C
- II. Re-suspend a single successful colony (determined from validating the transformation in 10mL liquid culture with antibiotic
- III. Incubate at 37°C until optical density reaches an absorbance of 0.4-0.8
- IV. Induce with 4 or 40µL of 100mM stock of IPTG (final concentration of 40 or 400µM) and induce for 3 to 5 hours at 37°C
- V. For large scale, inoculate 1L of liquid medium (with antibiotic) with freshly grown colony or 10 mL of freshly grown culture. Incubate at 37°C until optical density reaches an absorbance of 0.4-0.8. Add 40 or 400µM IPTG and express protein using optimal time and temperature determined in small scale trial

### Protein Purification under Denaturing

#### Ni-NTA Spin Kit

**NOTE:** due to dissociation of urea used pH values of buffers will need to be checked and if necessary adjusted

**NOTE:** this protocol is suitable for use with frozen cell pellets. Cell pellets frozen for at least 30 minutes at -20°C can be lysed by re-suspending in lysis buffer and adding Benzonase Nuclease (3units/mL culture volume). Fresh pellets require sonication or homogenization in addition. To the addition of 3 units/mL culture volume Benzonase Nuclease and 1mg/mL culture volume lysozyme.



- I. Thaw cells for 15 minutes and re-suspend in 700 $\mu$ L buffer B-7M urea and add 3 units/mL culture volume Benzonase Nuclease  
**NOTE:** Cells from 5 mL cultures are usually used, but culture volume used depends on protein expression level. Re-suspending pellet in 700 $\mu$ L buffer will allow recovery volume of cleared lysate of approx. 600 $\mu$ L
- II. Incubate cells with agitation for 15 minutes at room temperature. Solution should become translucent when lysis is complete.  
**NOTE:** buffer B is the preferred lysis buffer, as the cell lysate can be analyzed directly by SDS-PAGE. If the cells or the protein do not solubilize buffer A must be used.
- III. Centrifuge lysate at 12,000xg for 15-30 minutes at room temperature to pellet the cellular debris. Collect supernatant.  
**NOTE:** save 20 $\mu$ L of the cleared lysate for SDS-PAGE analysis
- IV. Equilibrate a Ni-NTA spin column with 600 $\mu$ L buffer B-7M urea. Centrifuge for 2 minutes at 890xg (2900RPM).  
**NOTE:** the spin columns should be centrifuged with an open lid to ensure that the centrifugation step is completed after 2 minutes.
- V. Load up to 600 $\mu$ L of the cleared lysate supernatant containing the GxHis-tagged protein onto a pre-equilibrated Ni-NTA spin column. Centrifuge for 5 minutes at 270xg (1600RPM), and collect the flow-through  
**NOTE:** to ensure sufficient binding, it is important not to exceed 270xg (1600RPM) when centrifuging Ni-NTA spin columns
- VI. Wash Ni-NTA spin column with 600 $\mu$ L buffer C. Centrifuge for 2 minutes at 890xg (2900RPM)
- VII. Repeat step VI
- VIII. Elute the protein with 200 $\mu$ L buffer E. Centrifuge for 2 minutes at 890xg (2900RPM) and collect the elute
- IX. Repeat step VIII

## SDS PAGE and Silver Staining

### Make PA Gel Separating Gel

- I. Level glass plates in the casting frame and place the frame within the casting stand. Ensure the plates are locked into place.
- II. Place the following solutions in a 50 mL falcon tube:

Mix:

- 4.2 mL dH<sub>2</sub>O
- 2.5 mL TRIS 4X Buffer
- 3.3 mL Polyacrylamide Solution.
- III. Mix and invert six times.
- IV. Add 33.3 $\mu$ L of APS (Ammonium Persulfate) into the falcon tube.
- V. Mix and invert six times.
- VI. Add 6.7 $\mu$ L of TEMED to the falcon tube.  
**NOTE:** Must be done in a fume hood. TEMED is toxic.
- VII. Mix and invert six times. Solution will harden in roughly 15 minutes.
- VIII. Add solution in between glass plates using a pipette aid roughly 3/4 the way up the front plate.
- IX. Fill remaining space with ethanol to ensure an even, level gel. Ethanol also gets rid of SDS bubbles.
- X. Once gel is solidified, remove the ethanol using a Kim Wipe ensuring the area in between the plates is now dry.

### Stacking Gel:

- I. Place the following solutions in a 50ml falcon tube:  
**Mix:**
  - 3.05 mL dH<sub>2</sub>O
  - 1.25 mL Stacking Buffer Solution
  - 650  $\mu$ L Polyacrylamide Solution
- II. Mix and invert six times.
- III. Add 25 $\mu$ L of APS (Ammonium Persulfate) into the falcon tube.
- IV. Mix and invert six times.
- V. Add 5 $\mu$ L of TEMED to the falcon tube.  
**NOTE:** Must be done in a fume hood. TEMED is toxic.
- VI. Mix and invert six times.
- VII. Fill in the remaining area in between the glass plates with the stacking solution and place comb inside.
- VIII. After gel is solidified, remove from casting frame and place in plastic wrap. Store in the fridge.
- IX. Gel is ready for use

### Run Electrophoresis

- I. Prepare samples that will be run in the gel electrophoresis.
- II. Mix 15 $\mu$ L of each sample with 5 $\mu$ L of sample buffer solution. A 3 to 1 ratio is used. A maximum of 30 $\mu$ L can be inserted into each well.
- III. Place glass plates in the electrode assembly and into the tank cell.
- IV. Fill tank with 1X electro buffer solution. Using a loading pipette and tips, load desired samples into wells along with a molecular weight ladder sample.
- V. Connect the lid to the tank and place the correct cables onto the tank lid.

Run gel at 120 volts for 60-90 minutes or until desired bands are about the length of the gel.

## Silver Stain Plus

- I. Once SDS-PAGE is complete, remove gels from the glass plates and remove the stacking gel layer using a paper towel.
- II. Place gel in a glass container, wash with 200ml distilled water for 5 minutes. Repeat the wash step 4X. (Shake gel back and forth in shaker) Decant water into waste container.
- III. Mix a solution with 90ml of 10% acetic acid and 90ml of 10% methanol in a clean beaker. Add 10ml of fixative enhancer solution into the beaker and pour over the gel. Let sit for 20 minutes on shaker, or up to a maximum of overnight, depending on the desired sensitivity of the gel.
- IV. Decant solution into a proper waste container.
- V. Wash gel once again with 200ml of distilled water for 5 minutes on shaker. Repeat 4X. Decant water into waste container.
- VI. Mix both parts of stain solution into two separate falcon tubes:

<b>Mix 1</b>	<b>Mix 2</b>
2.5 mL of Silver Staining Complex Solution	1.25g of Development Accelerator
2.5 mL of Reduction Moderator Solution	25mL of DH20
2.5 mL of Image Development Reagent	
17.5 mL of DH20	

- VII. Mix both solutions together until each are dissolved, pour over gel. Let sit for up to 20 minutes or until desired bands show.
- VIII. Decant solution into proper waste container.
- IX. Prepare a stop solution containing 90ml of 10% acetic acid and 90ml of 10% methanol. Pour over gel and shake for 15- 20 minutes.
- X. Decant solution into proper waste container.
- XI. Wash gel with 200ml of DH20 for 5 minutes on shaker. Repeat 2X.
- XII. Review results