

iGEM TU/e 2015

Biomedical Engineering

Eindhoven University of Technology

Room: Ceres 0.04

Den Dolech 2, 5612 AZ Eindhoven

The Netherlands

Tel. no. +31 50 247 55 59

2015.igem.org/Team:TU_Eindhoven

Protein Expression

Table of contents

Protein Expression	1	Protein Expression	3
	1.1	Materials	3
	1.2	Setup & Protocol	3

1 Protein Expression

Estimated bench time: 60 minutes

Estimated total time: 17 hours

Purpose: Protein expression of the bacteria.

It is essential to work near the Bunsen burner at all times.

1.1 Materials

- 1.5 ml cuvettes
- 2YT-medium
- Aluminum foil
- Antibiotic stock(s)
- Arabinose (20%)
- Cell Density Meter (OD600)
- Fresh culture of bacteria containing the right plasmid(s).
- Incubator
- IPTG (1 M)
- Non-natural amino acid (10 mM)
- Pipettes and tips
- Sterile culture tubes

1.2 Setup & Protocol

- Prepare a culture tube containing 7 ml 2YT and 7 µl of both antibiotic stocks (ampicillin and chloramphenicol).
- Transfer 100 µl of the small culture to this new tube.
- Grow the bacteria in the incubator at 37 °C and 250 rpm.
- After 90 minutes: measure the OD600.
OD measurement requires a blank measurement with 1 ml 2YT.
Pipette 1 ml of the culture in the cuvette and measure the OD600.
- Put the culture back in the incubator (37 °C and 250 rpm). Regarding the fact that a cell division cycle takes around 20 minutes, calculate the amount of time the culture needs to obtain an OD600 of 0.6. (the OD600 doubles after ±20 minutes)
- After the additional time: measure the OD600 again. Pipette 1 ml of the culture in the cuvette and measure the OD600.
- The amount left in the culture should be 5 ml. When the OD600=0.6 wrap the culture tube in aluminum foil (only if your protein is sensitive to light) and add
 - 5.62 µl IPTG (1 M)¹
 - 56.24 µl arabinose (20%)¹
 - 562.42 µl non-natural amino acid (10 mM)¹
 This makes the final concentration in the culture tube:
 - 1 mM IPTG
 - 0.2% arabinose
 - 1mM unnatural amino acid.
- Perform protein expression of ±15 hours at 18 °C and 250 rpm.

¹ It depends on your plasmid(s) what you need to add to initiate protein expression.