

Nano dpop
Nanoview use:

הפעלה של מכשיר הנתון למטה. הנה הוראות השימוש.

1. Remove lens cleaner/ dry kimwipe.
2. Turn on Nanoview.
3. Choose life science (1), DNA (1). Press OK.
4. Put 2 μ l of blank solution on the mark (middle of the cross) and press blank button (OA/100%T), repeat this step one more time. Between blanking, wash the lens with deionized water.
5. Clean lens with Kimwipes (both sides of lens).
6. Put 2 μ l of sample and press Sample button. (Clean lens between samples with deionized water).
7. You get a result in ng/ μ l.
260/280 value (DNA/protein) that should be 1.8-2.0.
260/230 value (DNA/salts) that should be 1.6-2.2.
8. Clean lens with deionized water and Kimwipes and return lens cleaner/dry kimwipe.
9. Turn off the device.
10. Cover nanoview.