


Experienced User Protocol

All spins at $\geq 12,000 \times g$, except as noted.



1 Harvest & lyse bacteria 5 minutes at 5000 rpm Bacterial culture

- Pellet cells from 1-5 ml overnight culture (1 ml from TB or 2xYT; 1-5 ml from LB medium). Discard supernatant.
- Resuspend cells in 200 μ l Resuspension Solution. Pipette up and down or vortex. and move to eppendorf.
- Add 200 μ l of Lysis Solution. Invert gently to mix. Do not vortex. Allow to clear for ≤ 5 minutes
- * Prior to first time use, be sure to add the RNase A to the Resuspension Solution.



2 Prepare cleared lysate

- Add 350 μ l of Neutralization Solution (S3). Invert 4-6 times to mix.
- Pellet debris 10 minutes at max speed. 10 min!



3 Prepare binding column

- Add 500 μ l Column Preparation Solution to binding column in a collection tube.
- Spin at $\geq 12,000 \times g$, 1 minute. Discard flow-through.



4 Bind plasmid DNA to column


- Transfer cleared lysate into binding column. (118 S3)
- Spin ~~30~~ 1 minute. Discard flow-through.



5 Wash to remove contaminants



אין צורך
Top 10
הוא
הוא

- (Optional (EndA+ strains only): Add 500 μ l Optional Wash Solution to column. Spin 30", 1 minute. Discard flow-through.)
- Add 750 μ l Wash Solution to column. Spin ~~30~~ 1 minute. Discard flow-through.
- Spin 1 minute to dry column.
- * Prior to first time use, be sure to add ethanol to the concentrated Wash Solution.



6 Elute purified plasmid DNA

- Transfer column to new collection tube. אין צורך
- Add 100 μ l Elution Solution. Spin 1 minute. \rightarrow Add 50 μ l HBW/UPW
- * If a more concentrated plasmid DNA prep is required, reduce the elution volume to a minimum of 50 μ l.



תלוי מה
הוא
הוא
הוא

after adding 50 μ l hot HBW, wait 2 min. and then:

- (1) according to noa:
 - centrifuge for 1 min at 2,000 xg
 - " " " " 17,000 xg

- (2) according to sara:
 - centrifuge for 2 min at 17,000 xg

Pure Plasmid DNA