

BABS UNSW iGEM Lab Protocol						
Procedure	Name		Lactococcus Transformation			
	Description		Transformation of Lactococcus cells using electroporation			
Document	Name	Mackenzie Labine-Romain		Date	7/07/15	Version 1
Requirements	Time					
	PPE		Gloves, Labcoat			
	Equipment		Electroporation cuvette (2mm gap) BioRad Gene pulser Ice			
	Materials		DNA sample dissolved in 10mM Tris-HCl-EDTA SGM17MC (SGM17 with 20mM MgCl ₂ and 2mM CaCl ₂) Selection antibiotic Streptococcal regeneration medium agar plates			
Step 1	Thaw competent cells on ice.					
Step 2	Mix 40µl of cells with 1µl of DNA dissolved in Tris-EDTA.					
Step 3	Transfer suspension to ice-cooled electroporation cuvette and expose to a single pulse (gene pulser set at 25µF and at 2.0kV).					
Step 4	Immediately after, mix suspensions with 0.96ml of ice-cold SGM17MC and leave on ice for 5 mins.					
Step 5	Dilute cells in SGM17MC and incubate at 30°C for 2h.					
Step 6	Spread 100µl portions on streptococcal regeneration medium plates containing 1µg/ml of antibiotic.					
Step 7	Enumerate transformants after 1-2 days of incubation at 30°C.					
Notes	Adapted from: Holo, H., & Nes, I. F. (1989). High-frequency transformation, by electroporation, of Lactococcus lactis subsp. cremoris grown with glycine in osmotically stabilized media. <i>Applied and Environmental Microbiology</i> , 55(12), 3119-3123.					
Version History						