

# PCR Purification Protocol

Modified from GeneJET PCR Purification Kit

## Preparation:

Dilute Wash Buffer (concentrated) by adding 5:1 volume of ethanol (96-100%).

## Procedure:

All centrifugations should be carried out in a tabletop microcentrifuge at  $>12000\times g$  (10 000-14 000 rpm).

1. Add a 1:1 volume of Binding Buffer to completed PCR mixture. Mix thoroughly. Check the color of the solution. A yellow color indicates an optimal pH for DNA binding. If the color of the solution is orange or violet, add 10  $\mu\text{L}$  of 3 M sodium acetate solutions (pH 5,2) and mix. The color of the mix will become yellow.
2. If the DNA fragment is  $\leq 500\text{bp}$ , add a 1:2 volume of 100% isopropanol. Mix thoroughly.
3. Transfer up to 800  $\mu\text{L}$  of the solution from step 1 (or optional step 2) to the GeneJET purification column. Centrifuge for 30-60 s. Discard the flow-through.
4. Add 700  $\mu\text{L}$  of wash Buffer (diluted with the ethanol) to the GeneJET purification column. Centrifuge for 30-60 s. Discard the flow-through and place the purification column back into the collection tube.
5. Centrifuge the empty GeneJET purification column for an additional 1 min to completely remove any residual wash buffer.
6. Transfer the GeneJET purification column to a clean 1.5 mL microcentrifuge tube. Add 50  $\mu\text{L}$  of water to the center of the GeneJET purification column membrane and centrifuge for 1 min.

Note. For low DNA amounts the elution volumes can be reduced to increase DNA concentration. Elution volumes less than 10  $\mu\text{L}$  are not recommended.

7. Discard the GeneJET purification column and store the purified DNA at  $-20^{\circ}\text{C}$ .