

iGEM TU/e 2015

Biomedical Engineering

Eindhoven University of Technology Room: Ceres 0.04 Den Dolech 2, 5612 AZ Eindhoven The Netherlands Tel. no. +31 50 247 55 59 2015.igem.org/Team:TU_Eindhoven

InterLab Study: PCR purification



Table of contents

	- /		
PCR purification	/ 1	PCR purification	3
	/ 1.1	Materials	3
	1.2	Setup & Protocol	3
	2	Nanodrop	3
	2.1	Materials	4
	22	Setup & Protocol	4

1 PCR purification

Estimated bench time: 30 minutes Estimated total time: 30 minutes

Purpose: Purifying the product obtained from a PCR reaction.

It is essential to work with gloves at all times to protect the DNA from DNase activity.

1.1 Materials

- Autoclaved Eppendorf tubes
- Autoclaved H₂O (nuclease free water)
- PCR tubes with PCR product
- · Pipettes and tips
- QIAquick PCR Purification Kit
- MiniSpin centrifuge

1.2 Setup & Protocol

- Mix 1 volume of the PCR sample with 5 volumes of buffer PB.
- Load the sample on a QIAquick spin column which is inserted in a collection tube (with a maximum of 800 µL per run).
- Centrifuge the sample for 1 minute at 13,400 rpm. Weight-balance the sample well.
- Discard the flow-through.
- Wash the sample with 750 µl of PE buffer and centrifuge for 1 minute at 13,400 rpm.
- Discard the flow-through.
- Dry spin the sample for 1 minute at 13,400 rpm.
- Transfer the spin column to a new autoclaved Eppendorf tube.
- Load 42 μl of autoclaved H₂O on the column (pipette drops in the middle of the membrane, do not touch the membrane). Incubate for 1 minute and centrifuge for 1 minute at 13,400 rpm.
- The resulting elution product will contain purified PCR product.

2 Nanodrop

Estimated bench time: -

Estimated total time: 5 minutes start-up, 2 minutes per sample

Purpose: Determine the concentration of DNA samples.

You are working with DNA, so it is essential to work with gloves at all times to protect your plasmids from DNase activity.

2.1 Materials

- Autoclaved H₂O
- DNA samples
- NanoDrop spectrophotometer
- Pipettes and tips

2.2 Setup & Protocol

- Start the NanoDrop spectrophotometer.
- Select the DNA measurement 'Nucleic Acid' in the NanoDrop menu.
- Clean the surface of the NanoDrop with dH₂O.
- Preform a calibration and blank measurement by entering one drop of 2 μ l autoclaved H_2O .
- Clean the surface again and place 2 µl per sample on the NanoDrop and measure the concentration. Write down the concentration (possibly on cryo-babies to stick on your DNA samples).