

## Prepare Gibson reaction

### 1) pMel-10 + pAQR1

	pMel-10 PCR	gRNA pAQR1	Mastermix	Sterile MQ water
Conc fragment (ng/ $\mu$ L)	31,9	2,5	X	X
Volume ( $\mu$ L)	3	4	10	3

### 2) Positive control

Prepare positive control: 10  $\mu$ L of control + 10  $\mu$ L of Gibson mastermix

## Gibson reaction

2 hours 50 degrees

### Transform to *E. coli* (5 tubes of *E. coli*)

- Transforms 2  $\mu$ L of Gibson assembly to *E. coli*: construct and positive control.
- Also transform 2  $\mu$ L of 4X diluted Gibson assembly: construct and positive control (4X dilution: mix 5  $\mu$ L of Gibson assembly mix with 15  $\mu$ L of sterile MQ water).
- Do a negative control: transform with 2  $\mu$ L of 6X diluted pMel-10 PCR (6X dilution: 1  $\mu$ L pMel-10 PCR + 5  $\mu$ L sterile MQ water). This will transform 10 ng of linear empty plasmid.
- Plate on total five LB + amp plates (2x construct + 2x positive control + 1x negative control)
- 37 degrees' overnight (max 16h)