

## iGEM2016 – Microbiology – BMB – SDU

<b>Project type:</b> Plastic <b>Project title:</b> Characterizing and optimizing PHB production <b>Sub project:</b> 1. Determining the intracellular amount of PHB	<b>Creation date:</b> 2016.09.26  <b>Written by:</b> Joel  <b>Performed by:</b> Joel Mario Vej-Nielsen & Jakob Rønning
---	--

### 1. SOPs in use.

SOP0001\_v1 – ON culture of *E. coli*

SOP0042\_v1 - Gas chromatography

### 2. Purpose.

Proteomics analysis of Staphylococcal pantothenate kinase II in *E. coli*.

### 3. Overview.

Day	SOPs	Experiments
1	SOP0001_v1	ON cultures of <i>E. coli</i> .
2	SOP0042_v1	Gas chromatography
	SOP0042_v1	Gas chromatography

## 4. Materials required.

### Materials in use

Name	Components (Concentrations)	Manufacturer / Cat. #	Room	Safety considerations
Benzoic acid		Sigma Aldrich	Chemical room	
PHB	pure	Sigma Aldrich		Handle in fume closet
Methanol		Sigma Aldrich	FKF	
$H_2SO_4$	100%		Chemical room	
Chloroform		Sigma Aldrich	Chemical room	
Distilled water				
LB media		Sigma Aldrich	Chemical room	

## 5. Other

### Primer to PCR dilution:

## 6. Experiment history.

Date (YY.MM.DD)	SOPs	Alterations to SOPs and remarks to experiments
16.10.10	SOP0001_v1 ON culture	50 ml of ON culture was prepared for two biological replicates of bacteria expressing pantothenate kinase II and two biological replicates of <i>E. coli</i> with pSB1C3 alone. Samples were grown for 72 hours, then lyophilized and kept on ice.
16.10.13	SOP0042_v1 Gas chromatography	The injection split ratio was 120:1. The injection port and the detector temperatures were set at 180°C and 200°C, respectively. The initial oven temperature was constant for one minute at 90°C.

		<p>Then increased at <math>3.8 \frac{^{\circ}\text{C}}{\text{minute}}</math> to a final temperature of 150°C.  The temperature was maintained for 5 min.  The flow rate of the helium carrier gas was 1 mL/minute.  First the analysis was performed with different concentrations of PHB and benzoic acid, in order to determine appropriate values.  Also Top 10 was analyzed as a negative control</p>
<b>16.10.13</b>	SOP0042_v1 Gas chromatography	<p>The injection split ratio was 120:1.  The injection port and the detector temperatures were set at 180°C and 200°C, respectively.  The initial oven temperature was constant for one minute at 90°C.  Then increased at <math>3.8 \frac{^{\circ}\text{C}}{\text{minute}}</math> to a final temperature of 150°C.  The temperature was maintained for 5 min.  The flow rate of the helium carrier gas was 1 mL/minute.</p>

## 7. Sample specification.

Sample name	Sample content	From	Used for / Saved where
#94	Top10 containing K2018036	Assembled from genscript genes and iGEM	Sample tested containing pantothenate kinase II
#9	Top10	Assembled from genscript genes and iGEM	Sample tested containing pantothenate kinase II

## 8. Remarks on setup.

## 9. Results and conclusions.

## 10. Appendixes