

Tar intein incubation protocol

Goal: To check the intein splicing

Protocol:

3 ml O/N starter in TB 30C to:

*Tar1 (+CM)

*UU

*Tar1 intein2 (+CM)

1. Dilute 1:50 with TB and growth in 30°C to O.D(600nm) ~0.6* (Erlenmeyer 250 ml) . **Total volume:** 35 ml. (700 ul Bactria+ 34.3 ml TB)
2. Dilution with the bacteria: 0, 10⁻¹, 10⁻², 10⁻³, 10⁻⁴, 10⁻⁵. (to falcon 50ml)
 0-w/o
 0-Etoh
 (7 tub per bacteria, 21 total)

	Cons. 10 ⁻¹	Cons. 10 ⁻²	Cons. 10 ⁻³	Cons. 10 ⁻⁴	Cons. 10 ⁻⁵
Bactria	4.95 ml	4.5 ml	4.5 ml	4.5 ml	4.5 ml
Bactria with β-estradiol	50 μl* β-estradiol fresh (C=10 mg/ml)	500 μl	500 μl	500μl	500μl

controls	with water	with absolute EtoH
bactria	4.95 ml	4.95 ml
water/EtoH	50 μl	50 μl

3. Move to 32C for 2 hr (for intein splicing).

-NOTE: Meanwhile make plats according to "swarming assay" protocol.

4. After 2 hr do "swarming assay".

*according to the article (the assay: "Chemical detection with strain DIER and the YES system").