

Bee T



Polymerase Chain Reaction

1. Keep everything on ice. The PCR reaction mixture in the ice preparation, and then placed in a PCR reaction in the PCR reaction apparatus. This cold start (Cool Start Method) can enhance the specificity of PCR amplification to reduce non-specific reactions in the PCR process, can get good PCR results.
2. Make up a master mix of everything into PCR tubes.
 - a. **25.0 μ L reaction system**

15.75 μ L ddH₂O
5 μ L 5X Q5 Reaction Buffer
0.5 μ L dNTPs
1.25 μ L forward primer
1.25 μ L reverse primer
0.25 μ L Q5 DNA polymerase
5 μ L 5X Q5 High GC Enhancer (optional)
1.0 μ L template DNA
-----**25.0 μ L** Total

b. 50.0 μ L reaction system

32.5 μ L ddH₂O
10 μ L 5X Q5 Reaction Buffer
1.0 μ L dNTPs
2.5 μ L forward primer
2.5 μ L reverse primer
0.5 μ L Q5 DNA polymerase
10 μ L 5X Q5 High GC Enhancer (optional)
1.0 μ L template DNA
-----**50.0 μ L** Total

3. Chose a suitable program, and adjust your annealing temperature and extention time as described below:

Annealing: 55°C for 0:30 min (different primers different annealing temperature)



Bee T



Extention: 72°C for t min (“t” depends on the length of goal sequence, 1min per 1000bp)

Final extension: 72°C for 10:00 min

