

iGEM2016 – Microbiology – BMB – SDU

Project type: Plastic Project title: Characterizing and optimizing PHB production Sub project: 1. Comparison of PHB producing strains	Creation date: 2016.08.24 Written by: Jakob Rønning Performed by: Joel Mario Vej-Nielsen & Jakob Rønning
---	---

1. SOPs in use.

SOP0001 – ON culture
SOP0042_v01 Daily use of FACS Aria
SOP0043_v01 Sorting setup FACS Aria II

2. Purpose.

Comparison of different PHB biobricks

3. Overview.

Day	SOPs	Experiments
1	SOP0001	ON culture
2	SOP0042	Sorting setup FACS Aria II
2	SOP0043	Daily use of FACS Aria

4. Materials required.

Materials in use

Name	Components (Concentrations)	Manufacturer / Cat. #	Room	Safety considerations
Appropriate medium ex. LB	1% Tryptone 1% NaCl 0.5% Yeast extract	Oxoid Sigma-Aldrich Merck	Media lab or V18-405-0	
Appropriate antibiotic if needed				
5 ml graduated pipettes		Fisher Scientific / CCI 4487	Micro storage	
Cuvettes		Contact lab-manager	BMB storage	
Blue pipette tips		Contact lab-manager	Micro storage	
Green pipette tips		Contact lab-manager	Micro storage	
Plate Bag		Contact lab-manager	BMB storage	
CST beads	Cytometer Setup and Tracking beads consist of bright (3 µm), mid (3 µm), and dim (2 µm) beads dyed with a mixture of fluorochromes that are excited by the lasers used in BD digital flow cytometers. See paragraph 12. Appendix	BD / 642412	Refrigerator V16-501b-2	
Accudrop beads	This product contains a single population of 6-µm particles. Every particle contains a fluorophore that is excited at 670 nm and emits at 750 nm. The particles are supplied in 1.5 mL of water containing 0.05% Tween® 20 and 2 mM sodium azide. See paragraph 12. Appendix	BD / 345249	Refrigerator V16-501b-2	

PBS	1x Phosphate Buffered Saline (PBS Buffer): 1 L distilled H ₂ O 8 g NaCl 0.2 g KCl 1.44 g Na ₂ HPO ₄ 0.24 g KH ₂ PO ₄	Sigma-Aldrich / 31434N J. T. Baker / 3040-01 MERCK / 1.06580.1000 MERCK / 1.04873.1000	V18-405-0 V18-405-0 V18-405-0 V18-405-0	Autoclave before use
FACS Flow sheet fluid	See paragraph 12. Appendix	BD / 342003	Refrigerator V16-501b-2	
FACS Clean	See paragraph 12. Appendix	BD / 340345	Refrigerator V16-501b-2	
FACS Rinse	See paragraph 12. Appendix	BS / 340346	Refrigerator V16-501b-2	
Deionized H₂O		V18-405-0 or Hallway storage (1. Floor)		
FACS tube		BD / 352054	V16-501b-2	
15 ml tube		Sarstedt	BMB storage	
FACS Flow sheet fluid	See paragraph 12. Appendix	BD / 342003	Refrigerator V16-501b-2	

5. Other

At the link <http://2015.igem.org/Team:Stanford-Brown/PHA>, "Using Nile Red to Quantify Plastic Production". We used their instructions to compare the different PHB biobricks in the FACS.

6. Experiment history.

Date (YY.MM.DD)	SOPs	Alterations to SOPs and remarks to experiments
-----------------	------	--



16.09.28	SOP0001	For #8, #9, #53, #90, #91, #92, #93 & #94 .
16.09.29	SOP0042	
16.09.29	SOP0043	

7. Sample specification.

Sample name	Sample content	From	Used for / Saved where
#8	<i>E. coli</i> with K2018000		Comparison
#9			Comparison
#53	<i>E. coli</i> with K934001		
#90	<i>E. coli</i> with K081005		
#91	<i>E. coli</i> with K2018033		
#92	<i>E. coli</i> with K2018034		
#93	<i>E. coli</i> with K2018035		
#94	<i>E. coli</i> with K2018036		

8. Remarks on setup.

9. Results and conclusions.



10. Appendixes

Data sheet for CST Beads: **23-9141-01_CS&T Beads_PI_RUO** is found in:

File name: iGEM2016_01_FACS_comparison_of_cells.docx

S:\999 Public\ISOPs and Protecols\FACS\FACS Aria II\SOPs for daily runs\Data sheets

Data sheet for Accudrop Beads: **23-8292-01_AccudropBeads** is found in:

S:\999 Public\ISOPs and Protecols\FACS\FACS Aria II\SOPs for daily runs\Data sheets

Data sheet for: FACS Flow sheet fluid: **342003-MSDS-EUEN-01 FACS Flow** is found in:

S:\999 Public\ISOPs and Protecols\FACS\FACS Aria II\SOPs for daily runs\Data sheets

Data sheet for: FACS Clean: **340345-MSDS-EUEN-00 FACS Clean** is found in:

S:\999 Public\ISOPs and Protecols\FACS\FACS Aria II\SOPs for daily runs\Data sheets

Data sheet for: FACS Rinse: **340346-MSDS-EUEN-00 FACS Rinse** is found in:

S:\999 Public\ISOPs and Protecols\FACS\FACS Aria II\SOPs for daily runs\Data sheets