



# **RISK ASSESSMENT – TASK BASED** **IGEM 2016**

<b>Location:</b> <i>Room W301, Medical Building</i>	<b>Building Number:</b> 181	<b>Date:</b> February 2016	<b>Assessed By:</b> Amber Willems Jones	<b>Health &amp; Safety Representative:</b> Vincé Kalangi
--	--------------------------------	-------------------------------	--	---

## **Description of Activity:**

### **4.8 Preparing Chemically Competent E.coli**

#### **SWP No: 4.8**

#### **Is there past experience with the Activity that may assist in the risk assessment?**

Incidents & Near-hits, Incident Investigations, Workplace Inspections, Training, Standards, Legislation & Codes, Uni Guidance Material, Existing Controls, Industry Standards.

**N  
O**

1. TASK	2. HAZARD	3. Estimated RAW RISK SCORE C x E x L	4. CONTROLS	5. Residual Risk SCORE C x E x L			6. Residual Risk	
Preparing Chemically Competent E.coli	Skin contact with Potassium chloride, Calcium Chloride, MOPS, MgCl2, liquid nitrogen, dry ice which may cause skin dryness and cracking, irritating, burns	25x3x1	Personal Protective Equipment ; training	25	3	0.1	7.5	low
Bench top centrifuge	Samples unbalanced, manual handling of rotors	15x3x1	Adequate training	15	3	0.1	4.5	
	TOTAL	120		TOTAL		12	low	
Name & Signature of Laboratory Head/Supervisor or Delegate	Amber Willems Jones		Date					
Name & Signature of Person Performing Activity or Task			Date					

<b>Number and Title</b>	PRG 4.8 Preparing Chemically Competent <i>E.coli</i>
<b>Name of Laboratory/Department</b>	The University of Melbourne IGEM Team Laboratory, Department of Biochemistry
<b>Author, Date Prepared &amp; Date of Review</b>	Author: Ella Bocquet-Gaylard Updated : February 2016, Date: 1/2/2016 Review by: February 2018
<b>Introduction</b>	The methods outlined in the following describe how to make competent cells
<b>Principles / Scope</b>	This method shows how to make E.coli cell competent for transformation
<b>Risk Management</b>	<b><i>Risk assessments have been prepared and are available attached to the SWP. Raw Risk: LOW Residual Risk: LOW</i></b>
<b>Safety Management</b>	<b>Hazards:</b> Always wear appropriate personal protective equipment. When handling hot materials. <b>Risk Controls:</b> Administrative, PPE
<b>Licences / Permits</b>	N/A
<b>Training / Competency</b>	All team members must be inducted into the use of any equipment used.
<b>Equipment</b>	Materials Sterile 50 ml Falcon Tubes Benchtop centrifuge Ice box Dry ice Liquid nitrogen Reagents TFB1 MOPS KCl CaCl <sub>2</sub> MnCl <sub>2</sub> Glycerol Hepes TFB2 HEPES KCl CaCl <sub>2</sub> Glycerol SOB medium Bacto-tryptone Yeast extract NaCl KCl MgCl <sub>2</sub> MgSO <sub>4</sub>

<b>Protocol</b>																																					
Step 1:	Grow 130 ml LB culture to OD <sub>600</sub> of 0.45.																																				
Step 2:	Pour cultures into sterile 50 mL falcon tubes.																																				
Step 3:	Pellet cells at 4 °C, 3000 rpm for 5 min.																																				
Step 4:	Discard supernatant.																																				
Step 5:	Resuspend pellet in total volume of 50 ml ice-cold TFB1.																																				
Step 6	Chill on ice for 5 mins.																																				
Step 7	Pellet cells at 4 °C, 3000 rpm for 5 mins.																																				
Step 8	Discard supernatant.																																				
Step 9	Resuspend pellet in 5 ml ice-cold TFB2.																																				
Step 10	Chill on ice for 15 mins.																																				
Step 11	Aliquot on dry ice, then snap-freeze, then store at -80 °C.																																				
	<b>Recipes</b> TFB1 (Transformation Buffer 1) components (500 mL): <table> <tr> <th>Reagent</th><th>Amount</th></tr> <tr> <td>10 mM MOPS</td><td>1.045 g</td></tr> <tr> <td>100 mM KCl</td><td>25 mL of 2 M sterile stock</td></tr> <tr> <td>10 mM CaCl<sub>2</sub></td><td>2.5 mL 2 M sterile stock</td></tr> <tr> <td>50 mM MnCl<sub>2</sub></td><td>12.5 mL of 2 M sterile stock</td></tr> <tr> <td>15 % Glycerol</td><td>75 mL</td></tr> </table> Dissolve MOPS and glycerol and pH to 5.8 with acetic acid and KOH in final volume of 460 mL then autoclave. Add the salts after autoclaving. TFB2 (Transformation Buffer 2) components (50 mL): <table> <tr> <th>Reagent</th><th>Amount</th></tr> <tr> <td>10 mM HEPES</td><td>0.11915 g</td></tr> <tr> <td>10 mM KCl</td><td>25 µL of 2 M sterile stock</td></tr> <tr> <td>75 mM CaCl<sub>2</sub></td><td>1.875 mL of 2 M sterile stock</td></tr> <tr> <td>15 % Glycerol</td><td>7.5 mL</td></tr> </table> Dissolve HEPES and glycerol and pH to 6.5 with KOH in final volume of 48.1 mL. SOB medium (500 mL) <table> <tr> <th>Reagent</th><th>Amount (g)</th></tr> <tr> <td>Bacto-tryptone</td><td>10</td></tr> <tr> <td>Yeast extract</td><td>2.5</td></tr> <tr> <td>10 mM NaCl (0.6 g/L)</td><td>0.3</td></tr> <tr> <td>2.5 mM KCl (0.18 g/L)</td><td>0.09</td></tr> <tr> <td>10 mM MgCl<sub>2</sub> (MgCl<sub>2</sub>·6H<sub>2</sub>O 2.03 g/L)</td><td>1.015</td></tr> <tr> <td>10 mM MgSO<sub>4</sub> (MgSO<sub>4</sub>·7H<sub>2</sub>O 2.46 g/L or anhydrous 1.35 g/L 1.35 g/L)</td><td>1.23/0.675</td></tr> </table> pH to 7.0 with NaOH, autoclave and store at 4 °C.	Reagent	Amount	10 mM MOPS	1.045 g	100 mM KCl	25 mL of 2 M sterile stock	10 mM CaCl <sub>2</sub>	2.5 mL 2 M sterile stock	50 mM MnCl <sub>2</sub>	12.5 mL of 2 M sterile stock	15 % Glycerol	75 mL	Reagent	Amount	10 mM HEPES	0.11915 g	10 mM KCl	25 µL of 2 M sterile stock	75 mM CaCl <sub>2</sub>	1.875 mL of 2 M sterile stock	15 % Glycerol	7.5 mL	Reagent	Amount (g)	Bacto-tryptone	10	Yeast extract	2.5	10 mM NaCl (0.6 g/L)	0.3	2.5 mM KCl (0.18 g/L)	0.09	10 mM MgCl <sub>2</sub> (MgCl <sub>2</sub> ·6H <sub>2</sub> O 2.03 g/L)	1.015	10 mM MgSO <sub>4</sub> (MgSO <sub>4</sub> ·7H <sub>2</sub> O 2.46 g/L or anhydrous 1.35 g/L 1.35 g/L)	1.23/0.675
Reagent	Amount																																				
10 mM MOPS	1.045 g																																				
100 mM KCl	25 mL of 2 M sterile stock																																				
10 mM CaCl <sub>2</sub>	2.5 mL 2 M sterile stock																																				
50 mM MnCl <sub>2</sub>	12.5 mL of 2 M sterile stock																																				
15 % Glycerol	75 mL																																				
Reagent	Amount																																				
10 mM HEPES	0.11915 g																																				
10 mM KCl	25 µL of 2 M sterile stock																																				
75 mM CaCl <sub>2</sub>	1.875 mL of 2 M sterile stock																																				
15 % Glycerol	7.5 mL																																				
Reagent	Amount (g)																																				
Bacto-tryptone	10																																				
Yeast extract	2.5																																				
10 mM NaCl (0.6 g/L)	0.3																																				
2.5 mM KCl (0.18 g/L)	0.09																																				
10 mM MgCl <sub>2</sub> (MgCl <sub>2</sub> ·6H <sub>2</sub> O 2.03 g/L)	1.015																																				
10 mM MgSO <sub>4</sub> (MgSO <sub>4</sub> ·7H <sub>2</sub> O 2.46 g/L or anhydrous 1.35 g/L 1.35 g/L)	1.23/0.675																																				
<b>Controls / Calibration</b>	N/A																																				
<b>Waste Disposal</b>	<b><u>Disposal requirements:</u></b> Follow PC I guidelines for handling, cleaning and when necessary, disposal of bacterial culture and solid wastes.																																				
<b>Emergency Procedures</b>	First aid measures Eye contact: Immediately flush eyes with plenty of water for at least 20																																				

	<p>minutes and get medical attention.</p> <p>Skin contact: In case of contact, immediately flush skin with plenty of water for at least 20 minutes.</p> <p>Inhalation: Move exposed person to fresh air. If not breathing, if breathing is irregular or if respiratory arrest occurs, provide artificial respiration or oxygen by trained personnel. Get medical attention.</p> <p>Ingestion: Wash out mouth with water. Do not induce vomiting unless directed to do so by medical personnel. Never give anything by mouth to an unconscious person. Call medical doctor or poison control centre immediately.</p>
<b>References</b>	
<b>Authorised By</b>	Amber Willems Jones