

Week 2

June 13, 2016

What We Did Today:

Made buffers

Split big flask into big 1:5 and into 2 small 1:20

Plated 12 well

Plating 12 Well:

First counts were too high so rediluted with 10 mL of media

New Counts:

237

322

Avg. 280

3×10^5 Calculations


$280 \times 4000 = 1.12 \times 10^6$ cells/mL

13 wells $\times 3 \times 10^5 = 3.9 \times 10^6$ cells in 13 mL

$3.9 / 1.12 = 3.48$ mL of cells

$13 - 3.48 = 9.52$ mL media added to cells

1 mL of media and cells was plated into each wells in the 12 wells plate



Handwritten calculations on a piece of paper:

$$225 \times 4000 = 900,000$$
$$3 \times 10^5 \times 13 = 3,900,000$$
$$\frac{3,900,000}{900,000} = 4.33 \text{ mL cell}$$
$$13 - 4.33 = 8.67 \text{ media}$$

Buffer creation:

2xHebs (200mL)

Mixed 2.00g of HEPES with 3.39g of NaCl

Dissolved in 190mL of sterile water

Adjusted pH to ~7.1 by adding drops of 1N NaOH

Adjusted volume to 200mL by adding sterile water

Filtered with .22 μ m disposable filter

70mM Na₂HPO₄ (200mL)

Dissolved 5.01g Na₂HPO₄·12H₂O with 200mL of sterile water
Filtered with .22µm disposable filter

2.5M CaCl₂ (200mL)

Dissolved 73.5g of CaCl₂·12H₂O in 200mL of sterile water
Stored at -20°C

1/10 TE buffer (1mM Tris-HCL and .1mM EDTA)

100 mL of sterile water
100 µL of 1M Tris-HCL
40 µL of .25M EDTA

General Lab Work

Transformation of pmCherry C3 and pTagRFP C1 into DH5α *E.coli* cells on KAN plates
Maintenance of HEK293T cell cultures

MG, JW

Split cell 1:4

Cell count:

294, 232 Average: 263

Plated in 12 wells at 3*10⁵ cells/mL

- Confluence: >100%, discard
- Split cell culture (Matthew's culture, 1:4)
- Cell count

Total count: count 1 = 294, count 2 = 232

- 4 mL cell + 9 mL media
- Plate in 12 wells at 3*10⁵ cells/mL

June 14, 2016

12 Well DNA Transfection Protocol

1.) CaPhos-DNA Mix

- a.) 50 μ L 2.5 M CaCl_2
 - b.) 10 μ g total DNA (5 μ g each plasmid)
 - I. I.) GFP: Use either pEGFPN3 \rightarrow 5 μ g = 7.7 μ L or pEGFPC1 \rightarrow 5 μ g = 5.7 μ L
 - II. II.) RFP: Use PTagRFPC \rightarrow 5 μ g = 5.7 μ L
 - c.) q.s. to 500 μ L with 0.1X TE
 - d.) Pipet up and down to mix
- 2.) Add 1 equal volume (500 μ L) of 2X HBS solution
 - 3.) Pipet up and down to mix (volume of tube is now 1 mL)
 - 4.) Incubate at room temp for 1 min
 - 5.) Add carefully 100 μ L each to 3 wells (this will be 1 μ g total DNA per well)
 - 6.) Add carefully 50 μ L each to 3 wells (this will be 0.5 μ g total DNA per well)

What We Did Today:

Transfected cells using GFP (pEGFPC1) and RFP

TouchTomorrow Data

Streaked InterLab colonies to new plates (Light Pink and Brown caps)

Sequenced rat APOBEC-XTEN-dCas9

Buffer Creation continued:

2xHBS (50mL) *all done in hood*

50mL of 2xHebs

1mL of Na_2HPO_4

Filtered with .22 μ m disposable filter

General Lab Work

Transfection of pTagRFP and eGFP into HEK293T cells

Maintenance of HEK293T cell cultures

Transformation of pmCherryC1 and pTagRFPC3 plasmids into 23716 *E.coli* cells

- Confluence: ~50%
- Transfection: pEGFP N3, pTag RFP
 - DNA 5 μ g each (7.7 μ L GFP, 5.7 μ L RFP)
 - Plated 100 μ L of mixture in 3 wells, 50 μ L in 3 wells
 - Change media after 4 hours

June 15, 2016

What We Did Today:

- Transformed pEGFPC1 and pEGFPN3
- Transformed InterLab Study brown and light pink cap (using 5 μ L DNA for each)
- Prepared new transfection plates
- Took pictures of Transfection plate with florescent microscope
- Made more art
- Split cells 1:5
- Named our team (RICE CRISPRS)

Cell Count Math

- Count = 250 cells
- $250 \times 4000 = 1,000,000$ cells/mL
- $3 \times 10^5 \times 13 = 3,900,000$
- $3,900,000 / 1,000,000 = 3.9$ mL of cells
- $13 - 3.9 = 9.1$ mL media

GM, JW

- Liquid culture: pmCherry C3, pTag RFP C1 (2 each)
- 5mL LB in each tube
- 5uL antibiotic (Kan) in each tube
- Pick colony and inoculate into conical tube
- Shake for a day

Made 100 mL cell culture media (JW)

- GM:
- Wells: 80% confluent
- 1:20 flasks: 40%
- 1:5 flash: frozen for stock

Cell counts:

- J: 228
- G: 241
- Avg: 234.5
- $\text{Avg} \times 4000 = 938000$
- $7 \times (3 \times 10^5) = 2.1 \times 10^6$ in 7mL
- $2.1 / .938 = 2.238$ mL of cells in 7mL
- $7 - 2.24 = 4.76$ media in 7mL
- 1mL per well 6 well

- Confluence: 90%
- Split culture
- Liquid culture: pmCherry C3, pTag RFP C1 (2 each)
- Cell count: 9.38×10^5 cells/mL
 - Plated in 6 wells at 3×10^5 cells/mL
- Made 100 mL new media

June 16, 2016


What We Did Today

- TouchTomorrow Report
- Poured Kan1 plates
- Transformed GFP plasmids
- Transfected cells with GFP and RFP
- Mini-Prep

Transfection for One Well

- 1.) 5 μL CaPhos
- 2.) 1 μg total DNA
- 3.) Dilute to 50 μL with 0.1x TE
- 4.) add 50 μL 2x HBS

Transfections for Today



GFP = 0.691 $\mu\text{g}/\mu\text{L}$

RFP = 0.87 $\mu\text{g}/\mu\text{L}$

M-Cherry = 0.54 $\mu\text{g}/\mu\text{L}$

Protocols Used:

75% GFP and 25% RFP

1. 1.) 15 μL CaPhos
2. 2.) 3.256 μL GFP

3. 3.) 0.862 μL RFP
4. 4.) Dilute to 150 μL with 0.1x TE
5. 5.) Add 150 μL 2x HBS

75% GFP and 25% m-Cherry

- 1.) 15 μL CaPhos
- 2.) 3.256 μL GFP
- 3.) 1.389 μL m-Cherry
- 4.) Dilute to 150 μL with 0.1x TE
- 5.) Add 150 μL 2x HBS

25% GFP and 75% RFP

- 1.) 15 μL CaPhos
- 2.) 1.085 μL GFP
- 3.) 2.586 μL RFP
- 4.) Dilute to 150 μL with 0.1x TE
- 5.) Add 150 μL 2x HBS

25% GFP and 75% m-Cherry

- 1.) 15 μL CaPhos
- 2.) 1.085 μL GFP
- 3.) 4.167 μL m-Cherry
- 4.) Dilute to 150 μL with 0.1x TE
- 5.) Add 150 μL HBS

GM, JW

Made glycerol stock of pTag C1 RFP and pmCherry

Mini-prep:

pmCherry: 17.6 ng/ μL
pTag: 32.3 ng/ μL (JW)
pmCherry 1: 54.0 ng/ μL
pmCherry 2: 87.5 ng/ μL

Transfection: 1 μg DNA/well

pmCherry only, 3 wells (GM)

pTag only, 3 wells (JW)

Liquid culture: pmCherry w/ kan added, 8 tubes

Transfections: (GM/JW)

Same protocol as before except:

Plated one row 100% pmCherry

one row 100% p-Tag

Transfection per one well:

5 μL CaPO₄

1 μg total DNA

Dilute to 50uL with .1xTE buffer
add 50uL 2xHBS

Use nanodrop to calculate total concentration and calculate off of that

pmCherry= .54ug/uL
1.85uL plasmid = 1ug DNA
pmCherry 3 wells:
15uL CaPO₄
5.56uL DNA
130uL .1xTE
150uL HBS

pTag = .87μg/μL
1.15μL plasmid = 1μg DNA
pTag 3 wells:
15μL CaPO₄
3.45μL DNA
131.55μL .1xTE
150uL HBS

- Glycerol stock of pTag C1 RFP
- Miniprep
 - pmCherry: 17.6 ng/μL
 - pTag: 32.3 ng/μL
 - (Discarded)
- Transfection: pTag RFP
 - 1 μg/well, plated in 3 wells
- Liquid culture:
 - pmCherry, 8 tubes with kan

June 17, 2016

What We Did Today

Mini-Prep m-Cherry plasmids

Split Cells

Photographed Transfections

Pet Dog

- Miniprep
pmCherry: 17.3 ng/ μ L
- Photo of 6/16 transfection: 50/50 w/ pmCherry, GFP, pTag, pmCherry
No significant bleed-over detected

GM:

Mini prepped pmCherry:

20ng/ μ L everyone got really low results (because of A4 buffer not containing ethanol)

Split 1:20 flasks to 1:10

forgot to bang one of the flasks after trypsin so some cells still stuck on flask

Took pictures of transfections

Checking for bleeding into other channels

Pictures on drive