

<h2>iGEM 2016 – SDU</h2>	
<b>Title:</b> Sample preparation	<b>Date issued:</b> 2016.10.04
<b>SOP number:</b> SOP0037	<b>Review date:</b> 2016.10.11
<b>Version number:</b> 01	<b>Written by:</b> Joel Vej-Nielsen

### 1. Purpose

Digesting cell pellet for phosphoprotein purification

### 2. Area of application

This procedure is valid for all *E. coli* strains.

### 3. Apparatus and equipment

Apparatus/equipment	Location (Room number)	Check points	Criteria for approval/rejection
Vacuum centrifuge			
Probe sonicator			
50ml falcon tubes			
Eppendorf tubes			
Centrifuge			
P20, P200 & P1000 pipettes			

### 4. Materials and reagents – their shelf life and risk labelling

Name	Components (Concentrations)	Manufacturer / Cat. #	Room	Safety considerations
PBS				
Lysis buffer (100µl)	93,8 µl Urea (6M) thiourea (2M) solution 1 µl 1M DTT (10mM) 3 µl phosSTOP phosphatase inhibitor 1,2 µl complete protease inhibitor 2 µl benzonase	Sigma Aldrich		Handle in fume closet
Triethylammonium bicarbonate (TEAB)	20mM, pH 7,5			
Iodoacetamide				
Trypsin				

## 5. QC – Quality Control

## 6. List of other SOPs relevant to this SOP

SOP0001 – ON culture  
SOP0038 – Qubit® Protein Assay Kits  
SOP0039 – C8 and C16 column purification  
SOP0040 - iTRAQ labelling  
SOP0041 - TiO2 purification

## 7. Environmental conditions required

## 8. Procedure

- 8.1 Prepare 50 ml cell culture of relevant bacterial strains.
- 8.2 Grow the cell cultures to an  $OD_{600} = 0.4-0.5$ .
- 8.3 Move 45 ml to a 50 ml falcon tube
- 8.4 Centrifuge the cell cultures for 5 min. at 4000 rpm and discard supernatant.
- 8.5 Resuspend cells in 10 ml PBS.
- 8.6 Centrifuge solution at 4000 rpm for 5 min. And discard the supernatant.
- 8.7 Resuspend cells in 5 ml PBS.
- 8.8 Centrifuge solution at 4000 rpm for 5 min. And discard the supernatant.
- 8.9 Resuspend cell pellet in 1ml PBS and move to Eppendorf tube.
- 8.10 Centrifuge the cell cultures for 5 min. at 4000 rpm and discard supernatant.
- 8.11 Store cells at  $-80^{\circ}\text{C}$ .
- 8.12 Resuspend appropriate amount of cell pellet in 15  $\mu\text{l}$  lysis buffer.
- 8.13 Incubate for 30 min. at room temperature.
- 8.14 Add 135  $\mu\text{l}$  TEAB.
- 8.15 Keep samples on ice and sonicate for 2x10 seconds.
- 8.16 Measure protein concentrations from 1  $\mu\text{l}$  of each sample on Qubit.
- 8.17 Add iodoacetamide to a final concentration of 20 mM and incubate for 30 min. in the dark.
- 8.18 Add 1  $\mu\text{g}$  of trypsin pr. 30  $\mu\text{g}$  protein and incubate over night at room temperature.
- 8.19 Dry down cells in speedy vac.

## 9. Waste handling

Chemical name	Concentration	Type of waste (C, Z...)	Remarks
Supernatant		Liquid bacterial waste	
Once use plastic		GMO yellow waste	

## 10. Time consumption

- Total-time 9h.
- Hands-on-time 1h.

## 11. Scheme of development

Date / Initials	Version No.	Description of changes
16.10.04 / JVN	01	The SOP has been written
16.10.11 / JR	01	The SOP has been reviewed

## 12. Appendixes