

## Risk assessment SB/IB-Rotary Shakers / Incubators

Produced 2010-05-11 By Riskbedömare 2 Livsvetenskaper at Systembiologi och industriell bioteknik Modified 2015-11-09 By Cecilia Sundelin

Final risk assessment of the	e method					
1. Acceptable risk						
1. State the premises in wh taking place	ich the activit	y is				
2. Description of activity						
3. Products						
Product name	Concentration	Form	Quantity	Danger	Comments	
Ammonium sulfate (mass)						
4. Risk category						
5. Level of exposure						
6. Ventilation						
7. Biological material						
8. Comments on Biological I	material					
9. Risk codes						
10. Comments on risk codes	5					
11. Premises						

12. Comments on premises

Risk assessment : SB/IB-Rotary Shakers / Inkubatorer

13. Protective signs
14. Comments on protective signs
15. Personal protective equipment
protective glasses, protective gloves, protective clothing
16. Comments on Personal protective equipment
Standard lab personal protective equipment
17. Describe the technical equipment
Rotary shakers are kept at a constant temperature (typically 30 or 37 oC) and are agitated to keep cells suspended and aid in aeration.
To use: 1 Stop shaker 2 Put shake flasks in proper position - On adhesive shakers, make sure bottoms are clean so they will stick well - On bar holders, tighten bar carefully, make sure they are straight. 3 Turn on shaker
18. Environment
19. Comments on environment
19. Comments on environment
20. Waste management
infectious radioactive or sharp waste , biological waste
21. Comments on Waste management
Shake flasks should be sterilized after use. Refer to shake flask protocol.
Shakers often have spilled cells on them, so consider the shake flasks contaminated at all times.  Glass shake flasks can break. Take care to clean up all the broken glass and be careful for sharp pieces of glass.  Boxes for contaminated glass waste in Analytical lab and Small lab, for clean broken glass in autoclave room.
22. Emergency equipment
fire-extinguisher foam, fire-extinguisher carbonic acid
23. Comments on Emergency equipment
Possible electric fire risk
Possible sharp, broken glass risk
24. Hazardous actions
night work

Risk assessment : SB/IB-Rotary Shakers / Inkubatorer

25. Comments on Hazardous actions
Do not be alone in the lab while operating orbital shakers.
26. Special instructions to other personel
27. Accidental readiness
28. Final risk assessment of the method
1. Acceptable risk
29. Comments on final risk assessment and additional risk reducing measurements
Signature
Supervisor Date
Christer Larsson

Date of reassessment: