

## Protocol for colony PCR

1.prepare the primer working solution, put Forward primer 1 $\mu$ l and Reverse primer 1 $\mu$ l into 8  $\mu$ l ddH<sub>2</sub>O

2.One kind of system(10 $\mu$ l):

Taq master mix: 5 $\mu$ l

Primer working solution: 0.1 $\mu$ l

DNA template: 0.1 $\mu$ l

ddH<sub>2</sub>O: 4.8 $\mu$ l

2.Another kind of system(10 $\mu$ l):

Hifi enzyme: 0.1 $\mu$ l

Primer working solution: 0.1 $\mu$ l

DNA template: 0.1 $\mu$ l

Hifi buffer: 1 $\mu$ l

dNTP: 1 $\mu$ l

ddH<sub>2</sub>O: 7.7 $\mu$ l

3.choose one kind of the two system above and mix them all together.

After preparation of the reaction mix, pick a single colony and dip in the mix. The reaction condition settings are as follows:

94°C 30s

55°C 30s

72°C 1kb/min

Cycles: 30-35