

COLONY PCR pfuX7

Sample preparation:

In the PCR tubes:

- 10µl MilliQ Water
- a small amount of the colony from the plate

After mixing the small amount of colony in the water, close the tubes and put them in liquid nitrogen for a few seconds. Immediately transfer them to a 37°C water bath to make a heat shock for 5 minutes. If no liquid nitrogen is available, put in PCR machine at 98°C for 10 minutes.

Preparation of the Master Mix:

PCR mix (on ice!)	1 x PCR mix (µl)
5x GC Buffer with MgCl ₂	4
2mM PCR Nucleotides	2
10µM Primer forward	1
10µM Primer reverse	1
2u/µl pfu X7	0.5
DMSO	1
MgCl	0.5

Add 10µl of the Mix in each PCR tube and place tubes in the PCR machine. Launch the appropriate program.

98° C	1'	<div style="display: inline-block; vertical-align: middle;"> <div style="border-left: 1px solid black; border-right: 1px solid black; height: 40px; width: 20px; margin: 0 auto;"></div> <div style="position: absolute; left: -5px; top: 50%; transform: translateY(-50%);">←</div> </div>	35 Cycles
98° C	20''		
52° C	45''		
72° C	2'30'		
72° C	5'		
Hold 12 ° C			

Spin down and use only the supernatant.