

<b>iGEM 2016 – SDU</b>	
<b>Title:</b> $TiO_2$ purification	<b>Date issued:</b> 2016.10.04
<b>SOP number:</b> SOP0041	<b>Review date:</b> 2016.10.11
<b>Version number:</b> 01	<b>Written by:</b> Joel Vej-Nielsen

## 1. Purpose

Procedure for isolation and enrichment of phosphopeptides in samples.

## 2. Area of application

This procedure is valid for digested protein samples.

### 3. Apparatus and equipment

Apparatus/equipment	Location (Room number)	Check points	Criteria for approval/rejection
SpeedyVac			
Vortex			
Centrifuge scale			

### 4. Materials and reagents – their shelf life and risk labelling

Name	Components (Concentrations)	Manufacturer / Cat. #	Room	Safety considerations
Loading Buffer	80% Acetonitril, 5% TFA & 1M glycolytic acid (76 mg/ml)			
Washing Buffer 1	80% Acetonitril & 1% TFA			
Washing Buffer 2	10% Acetonitril & 0,1% TFA			
Elution Buffer	60µl ammonia solution in 940µl $H_2O$ (pH 11,3)			
TFA	0,1%			
TFA	10%			
TFA	100%			
Acetonitril	30%			
HCl	12M			
HCl	1M			
PNGase F				
Sialidase A				
Formic acid	100%			
Low binding eppendorftubes				
$TiO_2$ beads				

## 5. QC – Quality Control

### 6. List of other SOPs relevant to this SOP

SOP0037 - iTRAQ sample preparation

SOP0038 – Qubit® Protein Assay Kits

SOP0039 – C8 and C16 column purification

SOP0040 - iTRAQ labelling

## 7. Environmental conditions required

## 8. Procedure

- 8.1 Mix the  $TiO_2$  loading buffer in your protein sample. (Add water, TFA Acetonitrile and glycolic acid in that order).
- 8.2 Add 0.6 mg  $TiO_2$  beads pr. 100 µg initial peptide solution.
- 8.3 Incubate on shaker (highest setting) for 15 minutes at room temperature.
- 8.4 Centrifuge sample briefly and move supernatant to another low binding Eppendorf tube with half the amount of  $TiO_2$  beads used previously.
- 8.5 Incubate for sample for 10 min. On shaker (highest setting).
- 8.6 Centrifuge briefly and move supernatant to another low binding Eppendorf tube (Mark tube “Non-modified peptides”).
- 8.7 Pool beads from both incubations with 100 µl binding buffer and transfer to new low binding Eppendorf tube
- 8.8 Vortex for 10 sec. and centrifuge briefly, before collecting supernatant in “Non-modified peptides” tube.
- 8.9 Wash beads with 100µl washing buffer 1, mix for 10 sec, centrifuge and move supernatant to “Non-modified peptides” tube
- 8.10 Wash beads with 100µl washing buffer 2, mix for 10 sec, centrifuge and move supernatant a new Eppendorf tube. (Mark tube “ $TiO_2$  wash2”)
- 8.11 Dry the  $TiO_2$  beads in vacuum centrifuge for 5-10 min.
- 8.12 Prepare 150 µl elution buffer.
- 8.13 Elute phosphopeptides by adding 150µl elution buffer and leave on shaker for 15. Min.
- 8.14 Centrifuge solution for 1 min. (highest speed) pass the supernatant over a filter (C8 stage tip) to recover liquid in a low binding tube.
- 8.15 Wash beads with 30µl elution buffer and pool supernatant with the other supernatant that has been filtered.

- 8.16 Elute from C8 filter with 5  $\mu$ l 30% acetonitrile
- 8.17 Add 7,3  $\mu$ l HCl (12 M) and adjust the pH of the solution to 7,5 with 1 M HCl solution.
- 8.18 Add 1  $\mu$ l PNGase F & 0.5  $\mu$ l sialidase A to the solution and incubate for one hour at 37°C.
- 8.19 Add 10  $\mu$ l 100% formic acid and 5  $\mu$ l 10% TFA to acidify solution.
- 8.20 Purify peptides by using two p200 R3 micro columns sequentially.
- 8.21 Wash columns with 60 $\mu$ l 0,1% TFA and elude phosphoproteins with 60  $\mu$ l 60% acetonitrile, 0,1% TFA. Start with the second column then the first one and then the eluate is collected in a new Eppendorf tubed (mark tube " $TiO_2$ -R3-eluate").
- 8.22 Run a MALDion 0,5  $\mu$ l sample and a LC-MS/MS on 5  $\mu$ l for individual group run.
- 8.23 Dry the samples by vacuum centrifugation
- 8.24 Lyophilize cells.

## 9. Waste handling

Chemical name	Concentration	Type of waste (C, Z...)	Remarks

## 10. Time consumption

- Total-time 5h
- Hands-on-time 3h

## 11. Scheme of development

Date / Initials	Version No.	Description of changes
16.10.04 / JVN	01	The SOP has been written
16.10.13 / JR		The SOP has been reviewed

## 12. Appendixes