

Yarrowia lipolytica LithAc transformation v1.0

Introduction

Transformation protocol from CAHOL, who works at CfB. Is optimized by .

Materials

- › MilliQ water
- › YPD plates
- › Appropriate selection reagents
- › YPD media
- › Micro centrifuge
- › 1.5mL eppendorf tubes
- › Either:
 - › Bürger-Türk hemocytometer with depth 0.1 mm squares of size C (0.2 mm x 0.2 mm)
 - › Microscope
- › Or:
 - › Cuvette
 - › OD meter
- › Transformation mix
 - › Polyethylene glycol (PEG)
 - › Lithium Acetate (LiAc)
 - › Salmon sperm DNA (ssDNA)
 - › Dithiothreitol (DTT)

Procedure

Prepare stock solutions for transformation mix

1. Prepare stock solutions using this protocol: http://openwetware.org/wiki/Springer_Lab:_TransformationYeast
2. Mix to prepare transformation mix

Make sure that all DTT is resuspended; it might be hard to see, but is crucial for good transformation efficiencies

Table1

	A	B	C	D
1		1x (μ L)	10x (μ L)	12x (μ L)
2	PEG (Stock 50%; sterile-filtrated; end 43.8%)	87.5	875	1050
3	LiAc (Stock 2M; sterile-filtrated; end 0.1 M)	5	50	60
4	ssDNA (Stock 10 mg/ml; end 0.25 g/l)	2.5	25	30
5	DTT (stock 2M; sterile-filtrated; end 100 mM)	5	50	60
6	Total	100	1000	1200

Day 1

3. Plate colony resuspended in 100 microliter MilliQ on YPD plates
4. Incubate plates for 24h @ 30°C

Day 2

5. Softly resuspended entire plate of cells in 1 mL MilliQ. Wash twice.
6. Centrifuge at 3000 RPM for 5-10 min.
7. Dilute cells 100-200x in MilliQ
8. Count using a Bürger-Türk hemocytometer with depth 0.1 mm. Count cells in 8 squares of size C (0.2 mm x 0.2 mm).
9. Calculate cells/mL = avarage of cells counted in 8 squares/volume of i square in mm³ * dilution factor x 10³
10. Otherwise: measure OD₆₀₀ and use a final OD of 9.2 per transformation
11. For each transformation use 5 x 10⁷ cells resuspended in 100 μ L transformation mix and 500ng DNA (max 5 μ L)
12. Incubate @39°C with shaking ~200RPM for 1h

01:00:00



13. If no antibiotics are used, plate the transformation mix directly on selective plates
14. If antibiotics are used, centrifuge the culture for 10 min at 3000 rpm and recovered in 600 μ L YPD for 2h before plating on selective plates.

02:00:00



15. Incubate @ 30°C for 2 days.

48:00:00

