iGEM 2016 – Microbiology – BMB – SDU

Project type: Bacteriocin

Project title: Cloning composite part into

iGEM standard plasmid.

Sub project: Insertion of k2018014

(Laterosporulin-Thuricin) into pSB1C3 plasmid

Creation date: 2016.08.09

Written by: Astrid Sophie Pejstrup Honoré

Performed by: Pernille Vigsø Rasmussen, Brian Kenn Baltzar, Astrid Sophie Pejstrup

Honoré, Cathrine Høyer Christensen

1. SOPs in use.

SOP number: SOP0007 v01 LA plates with antibiotic

SOP number: SOP0022 v01 Competent cell - freeze-stock

SOP number: SOP0023_v01 Ca⁺⁺ transformation

SOP number: SOP0009 v01 TSB transformation

Plasmid purification kit

SOP number: SOP0001_v01 ON culture of *E.coli*

SOP number: SOP0004 v01 bacterial freeze stock

SOP number: SOP0017 v01 Fast digest

SOP number: SOP0015 v01 Ligation

Gel purification kit

SOP number: SOP0021_v01 Colony PCR with MyTaq

2. Purpose.

To insert composite part: k2018014 (Laterosporulin-Thuricin), into iGEM standard plasmid; pSB1C3.

3. Overview.

Day	SOPs	Experiments
1	SOP0023_v01	Ca ⁺⁺ transformation
2	SOP0001_v01	ON culture of <i>E.coli</i>
3	Miniprep kit	Plasmid purification
3	SOP0004_v01	Bacterial freeze stock
4	SOP0017_v01	Fast digest
5	SOP0017_v01	Fast digest
5	Gel purification kit	Gel purification
6	SOP0015_v01	Ligation
7	SOP0009_v01	TSB Transformation
8	SOP0021_v01	Colony PCR with MyTaq
8	SOP0001_v01	ON culture of <i>E.coli</i>
9	Gel purification kit	Gel purification
9	SOP0004_v01	Bacterial freeze stock

4. Materials required.

Materials in use

Name	Components (Concentrations)	Manufacturer / Cat. #	Room	Safety considerations
Appropriate medium ex.	1% Tryptone 1% NaCl 0.5% Yeast extract	Oxoid Sigma-Aldrich Merck	Media lab or V18-40 5-0	
Glycerol	50 %	AppliChem	Anne Mette, RT	
LB		Anne-Mette		
LA	1% Tryptone 1% NaCl 0.5% Yeast extract 1.5% agar	Oxoid Sigma-Aldrich Merck Difco agar from BD	Anne-M ette Or	

V18-40 5-0

			5-0	
Water	Demineralised milli-Q autoclaved water	Milli-Q water purification system (Millipore)	RT	Water
MyTaq TM HS Red Mix	http://www.bioline.com/docu ments/product_inserts/MyTa q%E2%84%A2%20HS%20 Red%20Mix.pdf#zoom=130	Bioline	V18-405 a-2	
Reverse primer	Made specific to the template	Sigma-Aldrich		
Forward primer	Made specific to the template	Sigma-Aldrich		
Ligasebuffer		Agilent Technologies	Freezer at 1. Floor	
Ligase			Freezer 1. Floor	Ligase
FastDigest enzyme		Agilent Technologies	Freezer at 1. Floor	
Fast digest green / 10 x FastDigest Buffer		Agilent Technologies	Freezer at 1.	
CaCl ₂	0.1M		Chem room	
MgCl ₂	0.1M		Chem room	MgCl ₂
liquid nitrogen	liquid nitrogen	liquid nitrogen	liquid nitrogen	
Fast digest green		Agilent Technologies	Freezer at 1.	
6x DNA Loading Dye		GeneRuler	fridge floor 1	
Fort. LB		the new Anne-Mette	Autocla ve room	
Polyethylen e glycol (PEG) 3.350		Sigma Aldrich	Micro Chemic al room	

Dimethylsul foxid (DMSO)	Sigma Aldrich	Micro Chemic al room
Magnesiumc 1M hloride (MgCl2) 1M	The New Anne-Mette	Autocla ve

5. Other

As competent cells, LB and LA media was used by all parts of our project and not just this protocol the dates for use of these SOPs are not added. this comment deal with SOP number: SOP0007_v01 and SOP0022_v01

Gel Electrophoresis is set at 75 V for 30-45 minutes, dependent on the gel percent.

6. Experiment history.

Date (YY.MM.DD)	SOPs	Alterations to SOPs and remarks to experiments
16.07.14		resuspended the received 4µg of plasmid with 40µl dH2O
16.07.14	SOP0023_v01 Ca++ transformation	Their where used 2.5 µl instead of the normal 1µl
16.07.16	SOP0001_v01 ON culture of <i>E.</i> <i>coli</i>	
16.07.17	miniprep kit	50 μl were used for elution step
16.07.17	SOP0004_v01	Name: #31
16.07.19	SOP0017_v01 Fast digest	Fast digest enzyme(pst1) cut overnight
16.07.20	SOP0017_v01 fast digest	Fast digest enzymes(EcoR1) cut in 30 min.
16.07.20	Gel electrophoresis	1kb plus DNA ladder(Green) was used

16.07.20	Gel purification kit	30 μl were used for elution step
16.07.21	SOP0015_v01	At the ligation step following ratios of vector and DNA were
	ligation	used: 1:1, 1:5 and 1:10. For cPCR most cultures are taken from both ratio, and 3 replicates, to minimise the risk of getting religation.

Ratio	1:1	1:5	1:10
Ligase buffer	2µl	2µl	2µl
Ligase	1μ1	1μl	1µl
BG52	1µl	5µl	10µl
PG17	1µl	1μl	1µl
H ₂ O	16µl	11µl	6µl

16.07.22	SOP0009_v01 TSB transformation	5 μL plasmid is used for every sample

16.07.24 SOP0021_v01

cPCR with MyTaq

Segement	Step	Temperatu re	Duration
1	Initial denaturati on	98 °C	2 min
2	34 cycles	98 °C	10 sec
		72 °C	15 sec
		72 °C	15 sec
3	Final extension	72 °C	5 min
4	Keep the sample	12 °C	HOLD

	a a 1 d	
	cold	

16.07.24	SOP0001_v01 ON culture of <i>E.</i> <i>coli</i>	
16.07.25	Gel electrophoresis on cPCR product	1 kb DNA ladder(green) was used
16.07.25	miniprep kit	
16.07.25	Freeze stock	Name; #48

7. Sample specification.

Sample name	Sample content	From	Used for / Saved where
BR35	pUCIDT-AMP:k2018 014	IDT	transferred to <i>E. coli</i> / saved in coolbox
BR39	pUCIDT-AMP:k2018 014	IDT	Purifeid BR35, used to Cut with restriction enzymes and ligated with the vector/saved in coolbox and on Freeze stock (#31)
PR17	pSB1C3:k748002	iGEM	Vector cut with pst1 and EcoR1 / plastic cool box
BG52	k2018014	IDT	Cuted BR39with pst1 and EcoR / Saved in coolbox
BB55	pSB1C3:k2018014		vector and composite part ligated and transferred to <i>E.coli</i> . /Saved in coolbox
BR65	pSB1C3:k2018014		purified and transferred BB54 /

8. Remarks on setup.

9. Results and conclusions.

The Gel Photo below here shows the results from cutting the BR39 with Pst1 overnight and EcoR1 at 30 minutes, and after that called BG52 (At the photo named; 52). The composite part is 458 bp. It is possible from the Gel photo to assume that the DNA is cut correctly. And from this gel is the BG51 band with around 500 bp purified, for a ligation with PR17, the vector.



The Gel Photo under this, shows from the left the ladder; in this gel is used Green ladder. the following bands are BB55, 1, 2 and 4 there is the ligation of the vector PR17 and the insert BG52. From the photo it can be expect that all 3, samples in 55 are the vector with insert.



10. Appendixes