

| | |
|--|--------------------------------|
| iGEM2014 – Microbiology – BMB – SDU | |
| Title: Ligation | Date issued: 2013.06.19 |
| SOP number: SOP0015_v01 | Review date: 2015.09.13 |
| Version number: 01 | Written by: ASF |

1. Purpose

To ligate pieces of DNA

2. Area of application

Cloning

3. Apparatus and equipment

| Apparatus/equipment | Location (Room number) | Check points | Criteria for approval/rejection |
|---------------------|------------------------|--------------|---------------------------------|
| Vortex | Laboratory 1. Floor | • | |
| Pipettes (p20, p10) | Micro Storage | • | |
| | | • | |
| | | • | |
| | | • | |
| | | • | |
| | | • | |

4. Materials and reagents – their shelf life and risk labelling

| Name | Components | Supplier / Cat. # | Room (hallway storage) | Safety considerations |
|---------------------|------------|----------------------|------------------------|-----------------------|
| Purple pipette tips | | Contact lab-manager | Micro storage | |
| Green pipette tips | | Contact lab-manager | Micro storage | |
| Eppendorftubes | | Contact lab-manager | Micro storage | |
| Distilled water | | Contact lab-manager | Micro storage | |
| Ligasebuffer | | Agilent Technologies | Freezer at 1. Floor | |
| Ligase | | | Freezer 1. Floor | |

| | | | | |
|-------------|--|--|----------------------|--|
| DNA piece 1 | | | Refrigiator 1. Floor | |
| DNA piece 2 | | | Refrigiator 1. Floor | |
| | | | | |
| | | | | |
| | | | | |

5. QC – Quality Control

6. List of other SOPs relevant to this SOP

7. Environmental conditions required

8. Procedure

1. Prepare the ligation mixture and mix by pipetting up and down
2. Leave the mixture overnight at 16°C
- 2a. If there is no time leave the ligation solution at 22.5°C for 30 mins. Then denature the ligase at 65°C for 10min.
3. Use ligation solution for transformations

| Reagents | Volume |
|--|--|
| 10x T4 DNA ligase buffer | 2 µL |
| T4 DNA ligase (add last!) | 1 µL |
| PCR product (cut) of each brick which is to be ligated – or 1 part plasmid and 5 part bricks | 5 µL or 10 fmol Plasmid, 0, 10 and 20 fmol PCR |
| H2O | to reach a total volume of 20µL |

9. Waste handling

| Chemical name | Concentration | Type of waste (C, Z...) | Remarks |
|-----------------|---------------|-------------------------|------------------|
| One use Plastic | | GMO | Yellow GMO Trash |

10. Time consumption

- 3 Hours
- 1 Hour + Ligation overnight

11. Scheme of development

| Date / Initials | Version No. | Description of changes |
|------------------------|--------------------|-------------------------------|
| 13.06.19 / ASF | 01 | The SOP has been written |
| 13.06.26 /PRA | 01 | The SOP has been approved |
| | | |
| | | |
| | | |

12. Appendixes