

iGEM 2016 – Microbiology – BMB – SDU

Project type: Bacteriocin	Creation date: 2016.08.10
Project title: Cloning composite part into iGEM standard plasmid.	Written by: Astrid Sophie Pejstrup Honoré & Brian Kenn Baltzar
Sub project: Insertion of k2018015 (Lacticin Q-Lacticin Z) into pSB1C3 plasmid	Performed by: Pernille Vigsø Rasmussen, Brian Kenn Baltzar, Astrid Sophie Pejstrup Honoré, Cathrine Høyer Christensen

1. SOPs in use.

SOP number: SOP0007_v01 LA plates with antibiotic

SOP number: SOP0022_v01 Competent cell - freeze-stock

SOP number: SOP0023_v01 Ca⁺⁺ transformation

SOP number: SOP0010_vo1 Phusion PCR

Plasmid purification kit

SOP number: SOP0001_v01 ON culture of *E.coli*

SOP number: SOP0004_v01 bacterial freeze stock

SOP number: SOP0017_v01 Fast digest

SOP number: SOP0015_v01 Ligation

Gel purification kit

SOP number: SOP0021_v01 Colony PCR with MyTaq

2. Purpose.

To insert composite part: k2018015 (Lacticin Q-Lacticin Z), into iGEM standard plasmid; pSB1C3.

3. Overview.

Day	SOPs	Experiments
1	SOP0023_v01	Ca ⁺⁺ transformation
2	SOP0001_v01	ON culture of <i>E.coli</i>
3	Miniprep kit	Plasmid purification
3	SOP0004_v01	Bacterial freeze stock
4	Phusion PCR	Phusion PCR
4	SOP0017_v01	Fast digest
4	SOP0015_v01	Ligation
5	SOP0023_v01	Ca ⁺⁺ transformation
6	SOP0021_v01	Colony PCR with MyTaq
6	SOP0001_v01	ON culture of <i>E.coli</i>
7	Gel purification kit	Gel purification not on Gel but PCR product
7	SOP0004_v01	Bacterial freeze stock

4. Materials required.

Materials in use

Name	Components (Concentrations)	Manufacturer / Cat. #	Room	Safety considerations
Appropriate medium ex. LB	1% Tryptone 1% NaCl 0.5% Yeast extract	Oxoid Sigma-Aldrich Merck	Media lab or V18-40 5-0	
Glycerol	50 %	AppliChem	Anne Mette, RT	
LB		Anne-Mette		
LA	1% Tryptone 1% NaCl 0.5% Yeast extract 1.5% agar	Oxoid Sigma-Aldrich Merck Difco agar from BD	Anne-M ette Or V18-40 5-0	

Water	Demineralised milli-Q autoclaved water	Milli-Q water purification system (Millipore)	RT	Water
MyTaq TM HS Red Mix	http://www.bioline.com/documents/product_inserts/MyTaq%E2%84%A2%20HS%20Red%20Mix.pdf#zoom=130	Bioline	V18-405 a-2	
Reverse primer	Made specific to the template	Sigma-Aldrich		
Forward primer	Made specific to the template	Sigma-Aldrich		
Ligasebuffer		Agilent Technologies	Freezer at 1. Floor	
Ligase			Freezer 1. Floor	Ligase
FastDigest enzyme		Agilent Technologies	Freezer at 1. Floor	
Fast digest green / 10 x FastDigest Buffer		Agilent Technologies	Freezer at 1.	
CaCl ₂	0.1M		Chem room	
MgCl ₂	0.1M		Chem room	MgCl ₂
liquid nitrogen	liquid nitrogen	liquid nitrogen	liquid nitrogen	
Fast digest green		Agilent Technologies	Freezer at 1.	
6x DNA Loading Dye		GeneRuler	fridge floor 1	
dH ₂ O		Contact Lab-manager	Micro Storage	
5X Phusion HF or GC Buffer		Agilent Technologies	Freezer at 1. Floor	
dNTPs (10 mM each dNTP)		Agilent Technologies	Freezer at 1. Floor	

Primer #1 (10µM)		Freezer at 1. Floor
Primer #2 (10µM)		Freezer at 1. Floor
Phusion® Hig h-Fidelity DNA Polymerase	New England Biolabs	Freezer at 1. Floor
PCR/gel purification kit		Micro storage

5. Other

As competent cells, LB and LA media was used by all parts of our project and not just this protocol the dates for use of these SOPs are not added. this comment deal with SOP number: SOP0007_v01 and SOP0022_v01

Gel Electrophoresis is set at 75 V for 30-45 minutes, dependent on the gel percent.

6. Experiment history.

Date (YY.MM.DD)	SOPs	Alterations to SOPs and remarks to experiments
16.07.26		resuspended the received 4µg of plasmid with 40µl dH2O
16.07.26	SOP0023_v01 Ca++ transformation	Their where used 2.5 µl instead of the normal 1µl
16.07.27	SOP0001_v01 ON culture of <i>E. coli</i>	
16.07.28	miniprep kit	50 µl were used for elution step, made in 3 samples and are after Nanodrop, mixed into one, and a new nanodrop concentration is stated.

16.07.28	SOP0004_v01 Freeze stock	Name; #52
----------	-----------------------------	-----------

16.07.29 Phusion PCR This is made to replace a defective promoter, Primer; k2018015_IDT

H ₂ O	32.5µl
5 x HFbuffer	10µl
dNTPs	1µl
Primer; k2018015_IDT	2.5µl
Suffix-Primer	2.5µl
Polymerase	0.5µl

Program:

STEP	Temp °C	Time (m)
1	95	02:00
2	95	00:10
3	72	00:16
4	72	00:15
5	72	10:00
6	20	HOLD

Step 2-4 was repeated 34 times.

16.07.29	Gel electrophoresis	1kb plus DNA ladder(Green) was used, in a 2% gel and 5µl sample
----------	---------------------	---

16.07.29 SOP0017_v01 Fast digest with XbaI and PstI - 15 hours

16.07.29 SOP0015_v01
Ligation At the ligation step following ratios of vector and DNA were used: 1:0(negative control), 1:5 and 1:10, BR16(positive control). For cPCR most cultures are taken from both ratio, and 3 replicates, to minimise the risk of getting religation.

Ratio	1:0	1:5	1:10
Ligase buffer	2µl	2µl	2µl

	Ligase	1µl	1µl	1µl
	BG78	0µl	5µl	10µl
	BG17	1µl	1µl	1µl
	H₂O	17µl	12µl	7µl

16.08.01 SOP0023_v01 Their where used 2.5 µl instead of the normal 1µl
Ca++
transformation

16.08.02 SOP0021_v01 Made on 23 colonies
cPCR with
MyTaq

Segement	Step	Temperat ure	Duration
1	Initial denaturati on	98 °C	2 min
2	34 cycles	98 °C	10 sec
		72 °C	15 sec
		72 °C	15 sec
3	Final extension	72 °C	5 min
4	Keep the sample cold	12 °C	HOLD

16.08.02 Gel 1 kb DNA ladder(green) was used
electrophoresis
on cPCR product

16.08.02 SOP0001_v01
ON culture of *E.*
coli

16.08.03 miniprep kit Name; BR82

16.08.03 Freeze stock Name; #62

7. Sample specification.


Sample name	Sample content	From	Used for / Saved where
BR72	pUCIDT-AMP:k2018015	IDT	transferred to <i>E. coli</i> / saved in coolbox
BR74	pUCIDT-AMP:k2018015	IDT	Purified BR72, used to Cut with restriction enzymes and ligated with the vector/saved in coolbox and on Freeze stock (#52)
BG17	pSB1C3:k748002	iGEM	Vector cut with pstI and XbaI / cool box
BY77	k2018015	IDT	PCR product / Saved in coolbox
BG78	k2018015		BY77 cut with pstI and XbaI. /Saved in coolbox
BB79	pSB1C3:k2018015		vector and composite part ligated and transferred to <i>E.coli</i> . /Saved in coolbox
BR82	pSB1C3:k2018015		purified and transferred BB79 / freeze stock in <i>E.coli</i> (#62)

8. Remarks on setup.


9. Results and conclusions.

From the Gel photo under this part, is the result from the pPCR shown. This is made to replace a defective promoter, with the correct one. To solve this, is the Primer; k2018015_IDT used. As it can be seen from the photo, is first the ladder, here the green one. and then sample BR74. The gene is _____ bp. and it is therefore possible from the gel to assume that the pPCR is succeed.

Note: There is loaded 5 µl of each sample.



The Gel Photo under this, shows on the first Gel from the left the ladder; in this gel is used Green ladder. the following bands are BB79, 1 to 14, on the other gel is the Ladder, and then BB79 sample 15 to 23. BB79 is the ligation of the vector BG17 and the insert BG78. From the photo it can be expected that all the samples, except; 2 and 6, are perfect ligated. and one of the colonies are therefore used to make ON for the plasmid purification and Freeze stock.



10. Appendixes