

# End of the Note Book

**16 September 2016**


## Colonies PCR

Mix (25  $\mu$ L total volume reaction):

- 12.5  $\mu$ L DreamTaq PCR Mastermix 2X (ThermoScientific)
- 2.5  $\mu$ L 10X DreamTaq Green Buffer (ThermoScientific)
- 0.25  $\mu$ L Forward Primer (0.5  $\mu$ M) iG001
- 0.25  $\mu$ L Reverse Primer (0.5  $\mu$ M) iG002
- 9.5  $\mu$ L H<sub>2</sub>O filtrated (qs 25  $\mu$ L)

+ clones

## Agarose gel 1 %



## Precultures

Positive clones selected from the agarose gel results were picked and used to inoculate 3 mL of LBC (25  $\mu$ g/mL chloramphenicol) medium.


## Observation

➔ All the precultures have grown

- PCR colony

Mix:

- Dream Taq MM 2X	10 $\mu$ L	280 $\mu$ L
- Fw	1.5 $\mu$ L	42 $\mu$ L
- Rv	1.5 $\mu$ L	42 $\mu$ L
- DNA	1 $\mu$ L	-
- H2O	6 $\mu$ L	168




Just to remind:

Position E5: Mix without clone

Position E6: Mix with pLAC into pSB1C3

Position E7: Mix with mRF-HisTag into Psb1C3

Migration map:



➤ Expected Bands

LDH\_V2= ~1550


PCT\_V2=~1950

PHAC\_V3=~2050


PHAC\_V4=~2000

➤ Cycle

95 °C	3 min
95 °C	30 sec
50 °C	30 sec
72 °C	1.15 min
72 °C	5 min
10 °C	$\infty$



Gel Plan migration:



**20 September 2016**

Precultures

3 mL LB medium in 15 mL Falcon tube+ 1  $\mu$ L *P. putida* KT2440 from -80 °C stock

Incubation 30 °C overnight with shaking.

- PCR Gibson of pct\_V2, phac\_V3, phac\_V4, pSB1B3 HQ, pSB1C36 mRFP-HisTag

Mix:

Vf= 50  $\mu$ L

- 5X Q5 Reaction Buffer
- 10 mM dNTP
- 10  $\mu$ M Fw
- 10  $\mu$ M Rv
- DNA
- Q5 (0.2 U/ $\mu$ L)
- H2O

10 $\mu$ L	70 $\mu$ L
1 $\mu$ L	7 $\mu$ L
2.5 $\mu$ L	- $\mu$ L
2.5 $\mu$ L	- $\mu$ L
1 $\mu$ L	-
0.5 $\mu$ L	3.5 $\mu$ L
32.5 $\mu$ L	227 $\mu$ L

X7

- Solubilization of primers

For PR\_IG075: Add 209  $\mu$ L H2O for [100  $\mu$ M]

For PR\_IG076: Add 198  $\mu\text{L}$  H<sub>2</sub>O for [100  $\mu\text{M}$ ]

For PR\_IG077: Add 242  $\mu\text{L}$  H<sub>2</sub>O for [100  $\mu\text{M}$ ]

For PR\_IG078: Add 236  $\mu\text{L}$  H<sub>2</sub>O for [100  $\mu\text{M}$ ]

For PR\_IG079: Add 257  $\mu\text{L}$  H<sub>2</sub>O for [100  $\mu\text{M}$ ]

For PR\_IG080: Add 300  $\mu\text{L}$  H<sub>2</sub>O for [100  $\mu\text{M}$ ]

For PR\_IG081: Add 235  $\mu\text{L}$  H<sub>2</sub>O for [100  $\mu\text{M}$ ]

For PR\_IG082: Add 229  $\mu\text{L}$  H<sub>2</sub>O for [100  $\mu\text{M}$ ]

For PR\_IG083: Add 238  $\mu\text{L}$  H<sub>2</sub>O for [100  $\mu\text{M}$ ]

For PR\_IG084: Add 266  $\mu\text{L}$  H<sub>2</sub>O for [100  $\mu\text{M}$ ]

Aliquot 1/10 of each primer

- Gel purification of PCT2, PhaCv3, PhaCv4
- Purification

## 21 September 2016

- Gel Extraction pSB1C3 + mRFP-Histag

Gel weight= 300 mg

- PCR clean-up pct2, phac\_v3, phac\_v4

Nanodrop:

Genes	Concentration (ng/ $\mu\text{L}$ )	Ratio 260/280
PCT_V2	152	1.82
PHAC_V3	130.7	1.84
PHAC_V4	132.9	1.83
pSB1C3 IHQ	108.3	1.80
pSB1C3 Biobrick	14.7	1.74

- Gibson
  - Ratio Calculation

➔ Pct2-phac3 into pSB1C3

pSB1C3 IGHQ	0.92 $\mu\text{L}$
H <sub>2</sub> O	6.75 $\mu\text{L}$

Mix:

PCT2	1.04 $\mu$ L
PHAC3	1.29 $\mu$ L

➔ Pct2-phac4 into pSB1C3 IGHQ

Mix:

pSB1C3 IGHQ	0.92 $\mu$ L
H2O	6.76 $\mu$ L
PCT2	1.04 $\mu$ L
PHAC4	1.27 $\mu$ L

➔ Pct2-phac3 into pSB1C3 Biobrick

Mix:

pSB1C3 Biobrick	6.80 $\mu$ L
H2O	0.87 $\mu$ L
PCT2	1.04 $\mu$ L
PHAC3	1.29 $\mu$ L

➔ Pct2-phac4 into pSB1C3 Biobrick

Mix:

pSB1C3 Biobrick	6.80 $\mu$ L
H2O	0.89 $\mu$ L
PCT2	1.04 $\mu$ L
PHAC4	1.27 $\mu$ L

Add 10  $\mu$ L of HIFI DNA Assembly MM in each.

And 10  $\mu$ L of hifi with 10  $\mu$ L of positive control

Keep on ice!

➤ Incubation: 50 °C, 15 minutes then on ice

- Chemically Competent Cells Transformation

➔ Add 2  $\mu$ L of the chilled assembled to 25  $\mu$ L competent cells

➔ Waiting 30 min in ice

➔ Heat shock at 42 °C for 30 sec

➔ Transfer the tube on ice for 2 min

➔ Add 950  $\mu$ L LB for each

➔ Incubation: 60 min, 37 °C, 250 rpm

➔ Spread 100  $\mu$ L on the plate

➔ Incubation 37 °C overnight

## 23 September 2016

- Gibson plate observation
  - ➔ DH5-alpha pct2-phac3-pSB1C3 iGHQ 14 clones
  - ➔ DH5-alpha pct2-phac4-pSB1C3 iGHQ 13 clones
  - ➔ DH5-alpha pct2-phac3-pSB1C3 Biobrick 1 clones
  - ➔ DH5-alpha pct2-phac4-pSB1C3 Biobrick 5 clones

- PCR Colony

### Mix:

- DreamTaq Master Mix 2X
- PR\_IG001
- PR\_IG002
- DNA
- H2O

10 µL	350 µL
1.5 µL	52.5 µL
1.5 µL	52.5 µL
1 µL	-
6 µL	210 µL

→  
x35

### Cycle:

95 °C	3 min
95 °C	30 sec
50 °C	30 sec
72 °C	1.30 min
72 °C	2 min
10 °C	∞

↻  
x25

- Result Gibson

N°2 pct2 phac4 pSB1C3 iGEM HQ

N°1 pct2 phac4

N°5 pct2 phac4 pSB1C3 Biobrick

- Preculture in LB+ Cam (10 µg/mL)
- Miniprep of the cultures
- Gel Verification of Gibson PCR
  - ➔ Don't work

### Growth experiments

Stock in 50 % glycerol stored at -80 °C

Experiments performed in triplicates on 96 wells plate with LB medium as blank.

The plate was incubated at 30 °C with shaking during 20 h after inoculating with 1/1000 dilution of preculture.

The OD at 600 nm was read every 20 s.

The experiments were repeated 3 consecutive days from fresh precultures started the day before.

The results were sent to Imperial College IGEM team for collaboration along with the results of carbon source tests performed in July.

### **27 September 2016**

- Sequencing  
➔ WatchBox 1887819 (27.09.16)

### **28th September 2016**

- Gibson  
➔ Calculation for Ratio 1:1

#### Mix DNA 1:

pct2	0.62 µL
Phac3	0.76 µL
H2O	8.62 µL

#### Mix DNA 2:

pct2	0.66 µL
Phac3	0.81 µL
H2O	8.63 µL

#### Mix DNA3:

pct2	0.62 µL
Phac4	0.75 µL
H2O	8.63 µL

Mix DNA 4:

pct2	0.66 $\mu$ L
Phac4	0.80 $\mu$ L
H2O	8.54 $\mu$ L

➔ Calculation for Ratio 1:2

Mix 5:

Diluted pct2	1.04 $\mu$ L
Diluted Phac3	1.29 $\mu$ L
H2O	6.75 $\mu$ L
IGHQ Vector	0.92 $\mu$ L

➔ Mix 6:

Diluted pct2	1.04 $\mu$ L
Diluted Phac4	1.27 $\mu$ L
H2O	6.76 $\mu$ L
IGHQ Vector	0.92 $\mu$ L

➔ Mix 7:

Diluted pct2	1.04 $\mu$ L
Diluted Phac3	1.29 $\mu$ L
H2O	4.27 $\mu$ L
Biobrick Vector	3.40 $\mu$ L

➔ Mix 8:

Diluted pct2	1.04 $\mu$ L
Diluted Phac4	1.27 $\mu$ L
H2O	4.29 $\mu$ L
Biobrick Vector	3.40 $\mu$ L

➔ Put 15 min to 50 °C and after at -20 °C

**29th September 2016**

- Transformation

➔ 2  $\mu$ L product + 25  $\mu$ L of Dh5-alpha cells+ 950  $\mu$ L SOC Medium

➔ 1 h Recovery



- Culture plate LB Cam (10 µg/mL)
- Preparing plate

$$C_i = 50 \text{ mg/mL} \quad v_i = \frac{10 \text{ µg} \times 25 \text{ mL}}{50000} = 5 \text{ µL}$$

$$C_f = 10 \text{ µg/mL}$$

$$V_f = 25 \text{ mL}$$

- Purification PCR of gBlock fragment
  - Nanodrop
- Nothing ☹

### 30th September 2016

- Concentration / SpeedVac
- Verification Gel of purifications
- PCR amplification gBlock

### 4 October 2016

#### Gibson assembly 2<sup>nd</sup> step

pSB1C3 IGEM Headquarters 110.3 ng/µL

pSB1C3 mRFP 14.1 ng/µL

*Mix Gibson from the 1<sup>st</sup> step giving the operon:*

- 1) PCT2+PhaC3 65.5 ng/µL
- 2) PCT2+PhaC4 60.8 ng/µL
- 3) PCT2+PhaC3 63.4 ng/µL
- 4) PCT2+PhaC4 65.7 ng/µL

*Gibson ratio 1:1*

Vector amplified from pSB1C3-mFRP	1)	Vector	3.40 µL (50 ng)
		Insert (mix 1)	1.21 µL
		H2O (qs 10 µL)	5.40 µL
	2)	Vector	3.40 µL (50 ng)
		Insert (mix 2)	1.30 µL
		H2O (qs 10 µL)	5.30 µL

	3)	Vector	3.40 µL (50 ng)
		Insert (mix 3)	1.25 µL
		H2O (qs 10 µL)	5.35 µL
	4)	Vector	3.40 µL (50 ng)
		Insert (mix 4)	1.20 µL
		H2O (qs 10 µL)	5.40 µL
Linear vector pSB1C3 provided by IGEM Headquarters	5)	Vector	0.91 µL (100 ng)
		Insert (mix 1)	2.42 µL
		H2O (qs 10 µL)	6.68 µL
	6)	Vector	0.91 µL (100 ng)
		Insert (mix 2)	2.60 µL
		H2O (qs 10 µL)	6.49 µL
	7)	Vector	0.91 µL (100 ng)
		Insert (mix 3)	2.50 µL
		H2O (qs 10 µL)	6.60 µL
	8)	Vector	0.91 µL (100 ng)
		Insert (mix 4)	2.41 µL
		H2O (qs 10 µL)	6.68 µL

Add 10 µL NEB HiFi DNA Assembly Mastermix 2X

Incubation 50 °C, 15 min

Cool on ice

### Transformation

- Thaw 50 µL chemo-competent DH5-alpha *E. coli*
- Add 2 µL DNA (mix Gibson)
- Incubation 30 min on ice
- Heat shock: 30 s at 42 °C
- Incubation 2 min on ice
- Add 950 µL NEB medium (#B9020S)
- Incubation 1 h at 37 °C with shaking
- Centrifugation 7000 rpm 1 min
- Resuspend cells with rest of supernatant after throwing away the most of it
- Spread all the content on plates LBC
- Incubation overnight 37 °C

**7 October 2016**

Solubilisation of IDT gBlocks in 100 µL of nuclease-free water and short centrifugation

Incubation 50 °C during 20 min

Vortexing and then, short centrifugation again

### -PCR of gBlock genes

Mix:

- 10  $\mu$ L Q5 reaction buffer 5X
- 1  $\mu$ L dNTP mix
- 2.5  $\mu$ L 10  $\mu$ M Forward Primer (iG063)
- 2.5  $\mu$ L 10  $\mu$ M Reverse Primer (iG064)
- 1  $\mu$ L DNA template
- 0.5  $\mu$ L Q5 DNA polymerase
- 32.5  $\mu$ L Nuclease-free water (qs 50  $\mu$ L)

Mix for 5 reactions

*PCR program:*

- 1) 98 °C 30 s
- 2) 98 °C 30 s
- 3) 55 °C 30 s
- 4) 72 °C 1 min  
35 repeats of the steps bloc from 2) to 4)
- 5) 72 °C 2 min
- 6) 10 °C infinite

### Colony PCR on Gibson transformation

- 12.5  $\mu$ L DreamTaq Green PCR Master Mix 2X (Thermo Scientific, #K1081)
- 1.25  $\mu$ L 10  $\mu$ M Forward Primer (iG001)
- 1.25  $\mu$ L 10  $\mu$ M Reverse Primer (iG002)
- 10  $\mu$ L Nuclease-free water (qs 25  $\mu$ L)

+ clone

Mix for 23 reactions

*PCR program:*

- 1) 95 °C 5 min
- 2) 95 °C 30 s
- 3) 50 °C 30 s
- 4) 72 °C 3 min  
30 repeats of the steps bloc from 2) to 4)
- 5) 72 °C 10 min
- 6) 10 °C infinite

-Nanodrop of PCR genes purification of 3.09.16

Genes	Concentration (ng/ $\mu$ L)	Ratio 260/280
LDH_V2	138.2	1.88
PCT_V2	115.6	1.92
PHAC_V3	125.4	1.91
PHAC_V4	134.7	1.91

**8 October 2016**

➔ Problem with several PCR reaction tubes ➔ we do them again

Colony PCR on Gibson transformation

- 12.5  $\mu$ L DreamTaq Green PCR Master Mix 2X (Thermo Scientific, #K1081)
- 1.25  $\mu$ L 10  $\mu$ M Forward Primer (iG001)
- 1.25  $\mu$ L 10  $\mu$ M Reverse Primer (iG002)
- 10  $\mu$ L Nuclease-free water (qs 25  $\mu$ L)


+ clone

Mix for 9 reactions

*PCR program:*

- 1) 95 °C 5 min
  - 2) 95 °C 30 s
  - 3) 50 °C 30 s
  - 4) 72 °C 3 min
- 30 repeats of the steps bloc from 2) to 4)
- 5) 72 °C 10 min
  - 6) 10 °C infinite

## Agarose gel 1 %



## Digestion E+P

Mix: Vf= 100  $\mu$ L

DNA	50 $\mu$ L
Buffer 2.1	10 $\mu$ L
EcoR1	1.5 $\mu$ L
Pst1	1.5 $\mu$ L
H2O	37 $\mu$ L

## PCR clean-up

➤ Nanodrop

7 October		
Genes	Concentration (ng/ $\mu$ L)	Ratio 260/280
LDH_V2	37.7	2.04
PCT_V2	29.8	2.07
Phac_V3	35.7	2.04
Phac_V4	40.3	2.05

30 September		
Genes	Concentration (ng/μL)	Ratio 260/280
LDH_V2	8.8	2.49
PCT_V2	7.7	2.32
Phac_V3	11.0	2.68
Phac_V4	8.1	2.46
21 September		
Genes	Concentration (ng/μL)	Ratio 260/280
LDH_V2	7.9	2.54
PCT_V2	7.8	2.59
Phac_V3	7.0	2.35
Phac_V4	7.3	2.31

3 September		
Genes	Concentration (ng/μL)	Ratio 260/280
LDH_V2	16.3	2.20
PCT_V2	7.1	2.79
Phac_V3	8.2	2.30
Phac_V4	9.0	2.52

## 9 October 2016

Pouring of 46 LBC (10 μg/mL chloramphenicol) plates + 5 LBK (50 μg/mL kanamycin) plates + 2 LBS (50 μg/mL spectinomycin) plates

## 10 October 2016

Concentration measurements:

LDH 7.10: 37.7 ng/μL  
PCT 7.10: 29.8 ng/μL  
PhaC3 7.10: 35.7 ng/μL  
PhaC4 7.10: 40.3 ng/μL

LDH 3.09: 16.3 ng/μL  
PCT 3.09: 7.1 ng/μL  
PhaC3 3.09: 8.2 ng/μL  
PhaC4 3.09: 9.0 ng/μL

LDH 21.09: 7.9 ng/μL  
PCT 21.09: 7.8 ng/μL  
PhaC3 21.09: 7.0 ng/μL  
PhaC4 21.09: 7.3 ng/μL

LDH 30.08: 8.8 ng/μL

PCT 30.08: 7.7 ng/ $\mu$ L  
PhaC3 30.08: 11.0 ng/ $\mu$ L  
PhaC4 30.08: 8.1 ng/ $\mu$ L

### 11 October 2016

- Preparation 1L water+ 0.1 % HCOOH
- Test of solubility -> 10  $\mu$ L of standard in 90  $\mu$ L Solvent

#### ➤ Anthocyanin

- ➔ Isobutanol = not soluble
- ➔ DMSO = not soluble
- ➔ Ethyl Acetate = not soluble
- ➔ CH<sub>3</sub>CN 60 % = not soluble
- ➔ CH<sub>3</sub>CN 100 % = not soluble
- ➔ ETOH = not soluble
- ➔ HCL = ~soluble
- ➔ Methanol + Water + HCL = ~soluble
- ➔ Methanol+ HCL= ~soluble

#### ➤ Quercetin

- ➔ EtOH = soluble

#### ➤ Maldivin

- ➔ Soluble in Methanol HCL

Bettencourt Collaboration = we decided to do spectrometry masse instead of HPLC because it is more fast in time and the quercetin will be passed in negative mode but the maldivin and anthocyanin will be passed in negative mode.

### 12 October 2016

- Solution Stock
- ➔ At 2 mg/mL in EtOH 100 %

2 mg -> 1 mL EtOH

12.3 mg -> 6.15 mL EtOH

- Standard range in H<sub>2</sub>O mQ

- Quercetin

Isotopic mass = 302.042664 Da

Formula =  $C_{15}H_{10}O_7$

Solution Stock: [2000 mg -> 1000 mL]

- [20 mg/1000 mL]  
➔ 6  $\mu$ L in 594  $\mu$ L  $H_2O$  mQ -> Solution 1
- [10 mg/L]  
➔ 300  $\mu$ L of 1 + 300  $\mu$ L  $H_2O$  -> Solution 2
- [5 mg/L]  
➔ 300  $\mu$ L of 2 + 300  $\mu$ L  $H_2O$  -> Solution 3
- [1 mg/L]  
➔ 100  $\mu$ L of 3 + 400  $\mu$ L  $H_2O$  -> Solution 4
- [0.5 mg/L]  
➔ 250  $\mu$ L of 4 + 250  $\mu$ L  $H_2O$  -> Solution 5
- [0.1 mg/L]  
➔ 50  $\mu$ L of 5 + 200  $\mu$ L  $H_2O$  -> Solution 6

- Maldivin

Isotopic mass = 330.073955 Da

Formula =  $C_{17}H_{15}O_7$

Solution Stock: 50  $\mu$ L + 450  $\mu$ L MeOH-HCl 2 M -> [100 mg/L]

- [20 mg/1000 mL]  
➔ 100  $\mu$ L in 400  $\mu$ L  $H_2O$  mQ or  $CH_3CN$  60 % -> Solution 1
- [10 mg/L]  
➔ 250  $\mu$ L of 1 + 250  $\mu$ L  $H_2O$  mQ or  $CH_3CN$  60 %-> Solution 2
- [5 mg/L]  
➔ 250  $\mu$ L of 2 + 250  $\mu$ L  $H_2O$  mQ or  $CH_3CN$  60 %-> Solution 3



- [1 mg/L]  
➔ 100 µL of 3 + 400 µL H<sub>2</sub>O mQ or CH<sub>3</sub>CN 60 %-> -> Solution 4
- [0.5 mg/L]  
➔ 250 µL of 4 + 250 µL H<sub>2</sub>O mQ or CH<sub>3</sub>CN 60 %-> -> Solution 5
- [0.1 mg/L]  
➔ 50 µL of 5 + 200 µL H<sub>2</sub>O mQ or CH<sub>3</sub>CN 60 %-> -> Solution 6
- Plaque MS  
➔ 75 µL of each dilution + 75 µL Solution C

Solution C = 95 % acetate + 5 % carbonate ammonium 10 mM

- Vial  
➔ 200 µL of each dilution in insert into vial for test HPLC

### 13th October 2016

- Ligation 1:10 with different Vf= 25 µL or Vf= 20 µL

Length

- GB3 = 2038 bp
- IGHQ = 2037 bp

### DNA Mix:

7 October 2016

LDH + IGFQ (100 ng)	
Vector	0.92 µL
H <sub>2</sub> O	2.92 µL
LDH	16.2 µL

LDH + GB3 (50 ng)	
Vector	2.16 µL
H <sub>2</sub> O	9.77 µL
LDH	8.08 µL

PCT + IGHQ (75 ng)	
Vector	0.69 µL

H2O	0.78 $\mu$ L
PCT	20.09 $\mu$ L
PHAC3 + GB3 (75 ng) 1:5	
Vector	2.16 $\mu$ L
H2O	2.44 $\mu$ L
PHAC3	25.9 $\mu$ L

PCT + GB3(50 ng)	
Vector	2.16 $\mu$ L
H2O	4.46 $\mu$ L
PCT	13.39 $\mu$ L

PHAC3 + IGHQ (75 ng)	
Vector	0.69 $\mu$ L
H2O	1.45 $\mu$ L
PHAC3	17.85 $\mu$ L

PHAC3 + GB3(50 ng)	
Vector	2.16 $\mu$ L
H2O	5.95 $\mu$ L
PHAC3	11.90 $\mu$ L

PHAC4 + IGHQ (75 ng)	
Vector	0.92 $\mu$ L
H2O	2.05 $\mu$ L
PHAC4	21.12 $\mu$ L

3<sup>rd</sup> September 2016

LDH + GB3 (50 ng) 1:10	
Vector	2.16 $\mu$ L
H2O	0.89 $\mu$ L
LDH	18.68 $\mu$ L

21th September 2016

PHAC4 + GB3(50 ng)	
Vector	2.16 $\mu$ L
H2O	7.29 $\mu$ L
PHAC4	10.56 $\mu$ L

30th September 2016

PHAC3 + GB3 (75 ng) 1:5	
Vector	2.16 $\mu$ L
H2O	2.51 $\mu$ L
PHAC3	30.33 $\mu$ L

- Preparation of standard range for MS

20 mg/L; 10 mg/L; 5 mg/L; 1 mg/L

- Preparation of samples  
75 µL of supernatant of each + 75 µL of solution C.

#### 14 October 2016

- Transformation of the ligation product by heat shock  
3 µL of ligation product into DH5-alpha strain.  
Recovery 1 hour.  
spread on LB + chloramphenicol plates (25 µg/µL).

#### 15 October 2016

##### Digestion PhaC3\_only and PCT2\_only

- 30 µL IDT gBlock gene
- 5 µL buffer NEB CutSmart 10X
- 1 µL EcoRI
- 1 µL SpeI
- 13 µL nuclease-free water (qs 50 µL)

Incubation 1 h 37 °C

DNA clean & concentrator-5 kit (Zymo Research, D4003S)

DNA binding buffer added with 5:1 ratio

1 min incubation after adding 20 µL of nuclease-free water before centrifugation for elution

##### Ligation

	Quick ligase	T4 DNA ligase
pSB1C3 digested by E+S (from clone B → 13.3 ng/µL)	6 µL	6 µL
gene digested by E+S	3 µL	3 µL

<b>Buffer</b>	10 $\mu$ L Quick ligase buffer 2X (NEB)	2 $\mu$ L T4 DNA ligase buffer 10X (NEB)
<b>ligase</b>	1 $\mu$ L Quick ligase (NEB)	1 $\mu$ L T4 DNA ligase
<b>H2O (nuclease-free water)</b>	-	8 $\mu$ L

#### Transformation ligation mixes in pSB1C3

- Thaw 50  $\mu$ L chemo-competent DH5-alpha *E. coli*
- Add 2  $\mu$ L DNA (mix ligation)
- Incubation 20 min on ice
- Heat shock: 45 s at 42 °C
- Incubation 3 min on ice
- Add 450  $\mu$ L LB medium
- Incubation 1 h at 37 °C with shaking
- Centrifugation 7000 rpm 3 min
- Throw away supernatant
- Resuspend cells with rest of supernatant (around 50  $\mu$ L)
- Spread all the content on plates LBC (10  $\mu$ g/mL chloramphenicol)
- Incubation overnight 37 °C

#### Transformation ligation mixes in pSEVA 224, 424 or 2311

- Thaw 50  $\mu$ L chemo-competent DH5-alpha *E. coli*
- Add 2  $\mu$ L DNA (mix ligation)
- Incubation 20 min on ice
- Heat shock: 45 s at 42 °C
- Incubation 3 min on ice
- Add 450  $\mu$ L LB medium
- Incubation 1 h at 37 °C with shaking
- Centrifugation 7000 rpm 3 min
- Throw away supernatant
- Resuspend cells with rest of supernatant (around 50  $\mu$ L)
- Spread all the content on plates LBK (50  $\mu$ g/mL kanamycin) for pSEVA 224 or 2311  
or LBS (50  $\mu$ g/mL spectinomycin) for pSEVA 424
- Incubation overnight 37 °C

DNA PCR clean-up kit (NEB) for gBlock genes amplified by PCR

**16 October 2016**

- ➔ Our strain of DH5-alpha bacteria seems already resistant to spectinomycin ➔ test by preculture in 3 mL LBS (50 µg/mL spectinomycin) liquid medium

Concentration of amplified gBlocks after PCR clean-up

PhaC4: 86.6 ng/µL  
PhaC3: 120.6 ng/µL  
LDH2: 202.5 ng/µL  
PCT2: 154.6 ng/µL

Digestion by E+P of amplified genes

	PhaC3	PhaC4	PCT2	LDH2
DNA genes	7 µL	7 µL	5 µL	4 µL
Buffer NEB CutSmart	5 µL	5 µL	5 µL	5 µL
EcoRI	1 µL	1 µL	1 µL	1 µL
PstI	1 µL	1 µL	1 µL	1 µL
H2O nuclease-free	36 µL	36 µL	38 µL	39 µL

DNA clean & concentrator-5 kit (#D4003S Zymo Research)

Ratio 5:1 of DNA binding buffer

Incubation 1 min with 20 µL of nuclease-free water before centrifugation for elution

Precultures

2 mL LBC (25 µg/mL chloramphenicol) + clones from transformation of DH5-alpha *E. coli* with ligation mixes of pSB1C3 with PhaC3\_only or PCT2\_only

PCT\_only: 9 clones from ligation with T4 DNA ligase and 5 clones from ligation with Quick ligase

PhaC3\_only: 1 clone from Quick ligase and 2 clones from T4 DNA ligase

Incubation of the 17 precultures overnight at 37 °C with shaking.

**17 October 2016**

Miniprep

GenElute Plasmid Miniprep Kit (Sigma, PLN350-1KT)

3 clones for PhaC3 and 14 clones for PCT2 (see above)

Testing methods to detect the interest compound in HPLC

Samples run

- Preparation of samples
  - Vortex and samples centrifugation, 2 min, max speed.
  - Take 120 µL of supernatant of each and put them in each 1.5 tube.
  - Add 120 µL acetonitrile 100 %.
  - Centrifugation, 1 min, max speed.
  - Load 200 µL in the vials for HPLC.

**18 October 2016**

Colony PCR on the last trials of insertion of PhaC3, PhaC4, PCT2, LDH2, PhaC3\_only and PCT2\_only inside pSB1C3 and of insertion of PhaC4, PCT2 and LDH2 inside pSEVA224 and pSEVA2311 in process.

PCR Colonies:

20 µL mix per sample with:

- 10 µL Green Taq Master Mix 2X (#K1081 Thermo Scientific)
- 0.25 µL Forward primer iG001 (10 µM)
- 0.25 µL Forward primer iG002 (10 µM)
- 9.5 µL H<sub>2</sub>O (qs 20 µL)

+clone

*PCR program:*

- 1) 95 °C 5 min
- 2) 95 °C 30 s
- 3) 50 °C 30 s
- 4) 72 °C 2 min
- 30 repeats of the steps bloc from 2) to 4)
- 5) 72 °C 10 min
- 6) 10 °C infinite

HPLC Analysis (see Bettancout collaboration)

Sequencings of promoters isolated from 224 and 2311 plasmids and born by pSB1C3 came back positive.

**20 October 2016**

Gel electrophoresis with mixture from PCR on colonies obtained from PhaC3\_only and PCT2\_only ligations in pSB1C3.



\*the two samples chosen are indicated by an arrow

Results for the sequencing of the promoters from pSEVA224 and pSEVA2311 came positives.

**21 October 2016**

Samples (promoters from pSEVA224 and pSEVA2311 and genes PhaC3\_only and PCT2\_only) sent to iGEM registry and to sequencing for PhaC3\_only and PCT2\_only.

**24 October 2016**

Sequencing results for PhaC3\_oly and PCT2\_only came positives.