

SOB frozen competent cells (Yoram Riter lab) – yields ~50 220ul aliquots

Preparation:

- 1) 500 ml erlenmeyer for cell growth
 - 2) 89 ml DI H₂O autoclaved
 - 3) 200 ml SOB, autoclaved (need 120ml + a few ml for blanks):
 - 4gr bacto-tryptone
 - 1gr bacto yeast extract
 - 0.1gr NaCl
 - 0.038gr KCl
 - 4) 100 ml CCMB MADE SAME DAY, sterile filtered with 0.22um filter
(need 50ml):
 - 1.18gr CaCl₂ 2H₂O
 - 0.4gr MnCl₂ 4H₂O
 - 0.2gr MgCl₂ 6H₂O
 - 1ml KOAc 1M (pH 7.5) (on Orna's bench)
 - 10ml glycerol 100% autoclaved
- CHILL ON ICE.**
- 5) CHILL: glass pippettes, tips, 3 50ml falcons, eppendorfs or deep 96-well plate + cover,
metal holder for eppendorfs or tray for 96-well plate, sterile trough for dispensing
cells (OPTIONAL)
 - 6) fill liquid N₂

Protocol:

- 0) grow cells overnight in LB at 37C, 250rpm
- 1) set centrifuge temperature to 4C, get 3-4 buckets of ice
- 2) dilute 1:100 into 120ml SOB (+antibiotics) in 500ml flask, grow to OD₆₀₀=0.3 at 37C, 250rpm (this is a good time to prepare CCMB and chill it)
- 3) chill cells on ice for 10min (in falcons – after this NO FIRE)
- 4) pellet: spin 5000rpm, 10min, 4C, in 3 falcons, 40ml each
- 5) 1st resuspension:
 - on ice, 10ml ice-cold CCMB per falcon (~1/3 original volume)
 - using chilled glass pippettes. combine into 1 falcon, add another 10ml CCMB.
- 6) chill on ice for 20min
- 7) pellet: spin 5000rpm, 10min, 4C
- 8) 2nd resuspension:

- on ice, total of 10ml ice-cold CCMB
using chilled glass pippettes
- 9) chill 10min on ice
- 10) aliquot 220ul into precooled eppendorfs or precooled sterile
96-well plate with
cover using multipipettor, chilled tips and sterile trough (place tips in
freezer)
- 11) throw into liquid N₂ to freeze quickly
- 12) remove from liquid N₂ with latex-gloved hands quickly (can use lid
of -80 box),
move to -80C box and store at -80C
- use heat shock protocol to transform.**