

Gibson Assembly

Material

NEBuilder HiFi DNA Assembly Reaction Master Mix

Procedure

1. Set up the following reaction on ice:

Recommended Amount of Fragments Used for Assembly

	Recommended Amount of Fragments Used for Assembly		
	2-3 Fragment Assembly	4-6 Fragment Assembly	Positive Control
Recommended DNA Ratio	vector:insert = 1:2	vector:insert = 1:1	
Total Amount of Fragments	0.03-0.2 pmols X μ l	0.2-0.5 pmols X μ l	10 μ l
NEBuilder HiFi DNA Assembly Master Mix	10 μ l	10 μ l	10 μ l
Deionized H ₂ O	10-X μ l	10-X μ l	0
Total Volume	20 μ l	20 μ l	20 μ l

2. Incubate samples in a thermocycler at 50°C for 15 minutes (when 2 or 3 fragments are being assembled) or 60 minutes (when 4-6 fragments are being assembled). Following incubation, store samples on ice or at -20°C for subsequent transformation.
3. Transform competent E. coli cells with 2 μ l of the assembled product, following the transformation protocol.