

Table:

Expected product and size:	pA15A, 2 kb	
H <sub>2</sub> O (reaction volume - reagents added = )	34.5 µL	
Buffer: Phusion HF Buffer (5X)	10 µL	
Primer 1: 15_01	2 µL	
Primer 2: 15_02	2 µL	
DNA Template: pA15A miniprep, 1ng/µL dilution	1 µL	
Polymerase: Phusion	0.5 µL	
Total reaction volume:	50 µL	
PCR Program:	Program name: iGEM PHUSION	
	98°C	2 min
	98°C	30 s
	(anneal temp)°C	30 s
	72°C	(extend time) min
	repeat steps 2-4: 30 times	
	72°C	4 min
Additional info: <ul style="list-style-type: none"> <li>- (did you use GC buffer instead?)</li> <li>- (did you prepare the tubes and have someone else put them into the machine when it was free or when their PCRs were also ready to run?)</li> </ul>		

## Preparations:

- Fill ice bucket.
- Find and thaw primers and buffer.
- Make 1 ng/µL dilution of template if PCRing from a miniprep.

## Procedure:

- 1) Acquire and label caps of PCR tubes; place on ice.
- 2) Pipette water, buffer, primers, and template. (Use an (n+1)X master mix if more than 2 tubes, where n is the number of tubes.)
- 3) Start PCR machine: Temperature Measurement Mode > Algorithmic, Sample Volume > 50 µL, Hot Start? > YES, Hot Start Temperature > 98°C. Machine will preheat to 98°C and pause.
- 4) Add 0.5 µL Phusion polymerase per tube.

5) Replace tubes in machine and hit resume button to start thermocycling.

Assembled 2015 by Olya Spassibojko. Adapted from 2012. Source: iGEM 2015.