

Protocols

Yeast genomic DNA extraction with LiOAc-SDS

PROTOCOL FOR:

Extraction of genomic DNA from yeasts for PCR-based applications

Marko Lõoke, Kersti Kristjuhan, and Arnold Kristjuhan
Institute of Molecular and Cell Biology, University of Tartu, Tartu, Estonia

BioTechniques Protocol Guide 2012 (p. 395) doi 10.2144/000113785

Read the related article online: www.BioTechniques.com/article/113672

We have developed a quick and low-cost genomic DNA extraction protocol from yeast cells for PCR-based applications. This method does not require any enzymes, hazardous chemicals, or extreme temperatures, and is especially powerful for simultaneous analysis of a large number of samples. DNA can be efficiently extracted from different yeast species (*Kluyveromyces lactis*, *Hansenula polymorpha*, *Schizosaccharomyces pombe*, *Candida albicans*, *Pichia pastoris*, and *Saccharomyces cerevisiae*). The protocol involves lysis of yeast colonies or cells from liquid culture in a lithium acetate (LiOAc)–SDS solution and subsequent precipitation of DNA with ethanol. Approximately 100 nanograms of total genomic DNA can be extracted from 1×10^7 cells. DNA extracted by this method is suitable for a variety of PCR-based applications (including colony PCR, real-time qPCR, and DNA sequencing) for amplification of DNA fragments of ≤ 3500 bp.

Procedure

1. Pick one yeast colony from the plate or spin down 100–200 μ L liquid yeast culture ($OD_{600} = 0.4$). Suspend cells in 100 μ L 200 mM LiOAc, 1% SDS solution.
2. Incubate for 5 min at 70°C.
3. Add 300 μ L 96–100% ethanol, then vortex.
4. Spin down DNA and cell debris at $15,000 \times g$ for 3 min.
5. Wash pellet with 70% ethanol.

6. Dissolve pellet in 100 μ L H_2O or TE and spin down cell debris for 15 s at $15,000 \times g$.
7. Use 1 μ L supernatant for PCR.

Reagents

0.2 M lithium acetate (LiOAc) 1%
SDS solution
Ethanol, 96–100% and 70%.

Address correspondence to Arnold Kristjuhan,
Institute of Molecular and Cell Biology,
University of Tartu, Riia 23, Tartu 51010,
Estonia. e-mail: arnoldk@ut.ee

To purchase reprints of this article, contact:
biotechniques@fosterprinting.com

Disposable Cuvettes for DNA/RNA!

BRAND UV-Cuvette



- UV-transparent and disposable
- Use between 220-900nm
- 70 μ L minimum sample volume
- 10mm lightpath

**Free
Samples!**

- Replace quartz cuvettes for DNA, RNA and proteins
- Eliminate risk of cross-contamination from cuvette re-use
- Convenient & economical: no more cuvette washing!
- Simplify sample storage: round caps hold tight
- Available individually wrapped, DNA-, DNase-, & RNase-free



Lab Rats Trust BrandTech®!

BRANDTECH®
SCIENTIFIC, INC.

Toll Free (888) 522-2726
www.brandtech.com