

# Bee T



## Colony PCR

1. Keep everything on ice. The PCR reaction mixture in the ice preparation, and then placed in a PCR reaction in the PCR reaction apparatus. This cold start (Cool Start Method) can enhance the specificity of PCR amplification to reduce non-specific reactions in the PCR process, can get good PCR results.
2. Make up a master mix of everything into one microcentrifuge tube.
3. Pipette up and down in the microcentrifuge tube, drain 25 $\mu$ L or 50.0 $\mu$ L solution to each PCR tube.

### a. 25.0 $\mu$ L reaction system

18.375 $\mu$ L ddH<sub>2</sub>O  
5 $\mu$ L 5X OneTaq® Standard Reaction Buffer  
0.5 $\mu$ L dNTPs  
0.5 $\mu$ L forward primer  
0.5 $\mu$ L reverse primer  
0.125 $\mu$ L OneTaq DNA polymerase  
Colony stab (template DNA)  
-----25.0 $\mu$ L Total

### b. 50.0 $\mu$ L reaction system

36.75 $\mu$ L ddH<sub>2</sub>O  
10 $\mu$ L 5X OneTaq® Standard Reaction Buffer  
1.0 $\mu$ L dNTPs  
1.0 $\mu$ L forward primer  
1.0 $\mu$ L reverse primer  
0.25 $\mu$ L DNA polymerase  
Colony stab (template DNA)  
-----50.0 $\mu$ L Total

4. Pick colonies from plates, dilute in 10  $\mu$ L ddH<sub>2</sub>O and add 1  $\mu$ L into PCR tubes.
5. Run the "Colony PCR" program, and adjust your extension time as described below.

### The "Colony PCR" program


Initial denaturation: 95°C for 5:00min  
30 cycles of:  
94°C for 0:30 min  
55°C for 0:30 min (different primers different annealing temperature)  
72°C for t min ("t" depends on the length of goal sequence, 1min per 1000bp)



# Bee T



Final extension: 72°C for 10:00 min



iGEM Wageningen UR 2016. Send us an [email](#) or follow us on [Facebook](#)!  
For more information about iGEM, check out their [website](#) and [Twitter](#)!