

August 1 2017 : Characterization of M15 T7 lysate and M15 T7 lysate from cells grown on YPTG

1. Aim:

To investigate protein synthesis of M15 T7 cell lysates with and without Top10-GamS as well as compare the usual M15 T7 lysate with lysate made from M15 T7 cells grown on YPTG.

2. Materials:



- Lysate from Top10+GamS cells
- Lysate from M15 T7 cells
- Lysate from M15 T7 cells grown on YPTG
- Platerreader
- 384 well plate
- Linear GFP + T7
- Energy Solution
- Buffer A
- GamS
- Nuclease-free water




3. Procedure:

[DNA] initial		DNA quantity	Which lysate	Lysate quantity	Energy solution	Buffer A	GamS	H2O
Linear T7	DNA	0.336	M15-T7	2.5	2.5	2.5	0.3	1.86
Linear T7	DNA	0.336	M15-T7 + Top10	2.5	2.5	2.5	0	2.16
Linear T7	DNA	0.336	M15-T7 YPTG	2.5	2.5	2.5	0.3	1.86
Linear T7	DNA	0.336	M15-T7 YPTG + Top10	2.5	2.5	2.5	0	2.16
none		0	M15-T7	2.5	2.5	2.5	0	2.5
none		0	M15-T7 YPTG	2.5	2.5	2.5	0	2.5
none		0	Top10-GamS	2.5	2.5	2.5	0	2.5


All quantities in microliters. Set up platerreader to measure 150 repeats with two minutes between each repeat, at 37C.

4. Results:





Integrated chart with the average taken of the graphs before :



5. Conclusion:

The repeats are similar and enable us to draw the following conclusions:

The M15 T7 lysate works better with Top10-GamS than with GamS added in protein form (compare: it takes half the M15 T7 lysate, yet the average protein synthesis is higher. Top10-GamS does not have the T7 polymerase, so all the GFP synthesis is due the M15 lysate, whose productivity is enhanced).

There was an error with the M17 T7 YPTG lysate, maybe the wrong lysate was taken, as we see expression comparable to the negative control (nonexistent).

The M15 T7 YPTG part of this experiment should be repeated for conclusion.

Graphe: enlever les YPTG

Label:

M15-T7 + gamS

M15-T7 + Top10-gamS

Both control

