

Spectinmyocin PCR Protocol

Introduction

Materials

>

>

Procedure

Procedure + Materials

1 x 50 μ L reaction (the gel I gave you is for a 4 x 50 μ L reaction)

- added water first and polymerase last, the order of all other reagents is irrelevant
- thawed DNTP's and buffer on ice
- kept the phusion polymerase in the freezer until the reaction was fully prepared, and then kept it on ice until use
- mixed everything (except water) thoroughly by flicking tubes

31 μ L n.f. water
10 μ L Phusion high-fidelity (HF) buffer
2.5 μ L JTL035 (forward primer)
2.5 μ L JTL036 (reverse primer)
1.5 μ L DMSO (optional, but I use 3% v/v of DMSO in all of my PCRs because of my poor annealing regions in my project, can replace this with an extra 1.5 μ L of n.f. water)
1 μ L of P5trc (name of the plasmid, concentration is roughly 5ng/ μ L)
1 μ L 10mM DNTP's
0.5 μ L of Phusion polymerase

Total volume: 50 μ L

Thermo-cycler conditions

1. Initial denaturation: 98°C, 30s

(iterate the next step 35 times)

2a. Denaturation: 98°C, 5s

2b. Annealing: 62.9°C, 15s

2c. Elongation 72°C 30s

3. Final elongation: 72°C 10min

4. Hold: 4°C, forever (until you remove the samples, keep at 4°C)

Gel Electrophoresis

Used 1% agarose gel, poured 1x TAE over the gel to fully submerge it, and loaded 50 μ L into each well (if you aim to PCR purify, only load 5 μ L) along with a 1x 1kb ladder in the first well. Ran at 100V for 30 minutes (you might have to run it longer because I remember our machines were wonky). Image gel and if gel extracting, excise correctly sized bands. **Correctly sized band: ~1kb**