

Protocol for Bacteria Transformation

.....

Aim

To transform competent bacteria with plasmids.

Material

- 500 µl tube of competent bacteria (competent cells stored at -80°C)
- Plasmid to transform
- 1.5 ml Eppendorf tubes
- P1000, P200, P10 pipettes + paired cones
- Petri dishes with LB agar/CARB for carbenicillin (an equivalent of ampicillin, resistant to temperature)
- Water-bath at 42°C
- Incubator at 37°C

Protocol

Plasmid transformed

-
-
-
-

Competent cells are extremely sensitive; always handle the tubes on ice.

You must operate in the vicinity of the Bunsen burner when manipulating bacterial cultures.

Name the tubes with the transformed cells: Cell type/name, plasmid (vector and composition), initials of the operator: First Name/Last Name, date.



Protocol for Bacteria Transformation

.....

1. Split the tube of competent cells in aliquots of 50 μ l, in 1.5 Eppendorf tubes placed on ice.
For maximum competency:
2. Add 1 μ l of β -mercapto-ethanol only to subcloning grade competent cells DH5alpha.
3. Mix by gently tapping the bottom of the tubes.
4. Add μ l = pg of plasmid to transform.
5. Mix by gently tapping the bottom of the tubes.
6. Let it rest for 30 minutes on ice.

In the meantime, check that the water-bath is at 42°C, and place the SOC media at 37°C for warming.

7. Put the tubes in the floats.
8. Place the floats in the 42°C water-bath for 40 s, then remove the floats quickly.
9. Place the tubes on ice for 3 min.
10. Add 650 μ l of SOC media per tube.
11. Incubate and mix the tubes at 150 rpm at 37°C for 40 min.

In the meantime, place the labelled LB/AMP petri dishes at 37°C.

12. Generate 2 petri dishes with each tube: one dish containing 200 μ l and the other 500 μ l.
13. Spread the bacteria using an inoculator.
14. Wait for the dishes to dry.
15. Store the dishes with agar side up in the incubator at 37°C for approximately 16 h (overnight).