

SDS-Page

- ◆ For each gel (12.5 %) aliquot:
 - ◇ 0.8 mL Tris-HCL (1.88 M), pH = 8.8
 - ◇ 0.8 mL SDS (0.5 %)
 - ◇ 1.67 mL Acryl-/Bisacrylamid (5:1)
 - ◇ 0.73 mL H₂O
 - ◇ 2.7 µL TEMED (99 %)
 - ◇ 20 µL APS (10%)
- ◆ Pour the solution quickly into the gel casting form. Leave about 2 centimeters below the bottom of the comb for the stacking gel.
- ◆ Layer isopropanol on top of the gel.
- ◆ Leave the separating gel at room temperature for >60 minutes to polymerize.
- ◆ Remove isopropanol and wash surface of the separating gel with H₂O. Wait until the surface is dry.
- ◆ For each stacking gel (5 %) aliquot:
 - ◇ 235 µL Tris-HCL (0.625 M), pH = 6.8
 - ◇ 235 µL SDS (0.5 %)
 - ◇ 195 µL Acryl-/Bisacrylamid (5:1)
 - ◇ 335 µL H₂O
 - ◇ 1.2 µL TEMED (99 %)
 - ◇ 6 µL APS (10 %)

- ◆ Insert comb without getting bubbles stuck underneath.
- ◆ Leave the gel at room temperature for > 60 minutes to polymerize.
- ◆ For storage:
 - ◇ Remove sealing and store the gel wrapped in moistened paper towel at 4 °C.
- ◆ Preparing the sample:
 - ◇ Mix your protein mixture 3:1 with PBJR buffer (15 µL protein solution + 5 µL PBJR buffer).
 - ◇ Heat for 5 minutes at 95 °C.
- ◆ Running the gel:
 - ◇ Remove sealing, put the polymerized gel into gel box and pour SDS-PAGE running buffer into the negative and positive electrode chamber.
 - ◇ Remove comp without destroying the gel pockets.
 - ◇ Pipet the SDS running buffer in the gel pockets up and down for flushing the gel pockets.
 - ◇ Pipet slowly 20 µL of the sample into the gel pockets.
 - ◇ Make sure to include at least one lane with molecular weight standards to determinate the molecular weight of the sample.
 - ◇ Connect the power lead and run the stacking gel with 10 mA until the blue dye front enters the separating gel.
 - ◇ When the distance of the lowest molecular weight standard lane to the gel end is down to 0.5 cm stop the electrophoresis by turning off the power supply.
- ◆ Staining the polyacrylamide gel (Colloidal Coomassie Brilliant Blue staining):
 - ◇ After finishing the SDS-PAGE remove gel from gel casting form and transfer it into a box.
 - ◇ Add 100 mL of the Colloidal Coomassie Brilliant Blue staining solution to your polyacrylamide gel.
 - ◇ Incubate the gel in the solution at room temperature until the protein bands got an intensive blue color. Shake the gel continuously during incubation.
 - ◇ Remove the staining solution.

- ◇ Wash the gel with 7 % (v/v) acetic acid in H₂O for decoloration.
- ◇ Incubate the gel in H₂O (2-6 h) for bleaching the background. Shake the gel continuously during incubation. If necessary, replace the colored water with new one.

From [Bielefeld-CeBiTec 2014](#)