

IMPACT-CN

- ◆ Preparation of Plasmids
 - ◇ Clone your protein fragment into one of the following vectors:
 - pTXB1
 - pTYB1
 - pTYB2
 - pTYB11
 - pTYB12
 - ◇ transfer the complete construct into ER2566 via electro transformation
- ◆ cell growth
 - ◇ prepare a pre-culture with some cell material within up to 10mL LB substituted with 200 µg/mL Amp and 20 µg/mL Cm
 - ◇ incubate the cultures overnight at 37°C
 - ◇ prepare the main culture:
 - 250mL LB (+Amp and Cm)
 - Inoculate the culture to an OD₆₀₀ of 0.1 with the pre culture
 - ◇ Incubate the main culture at 37°C until an OD₆₀₀ auf 0.5 to 0.8 is reached
 - ◇ Induce the protein expression via addition of IPTG (0.3-0.5mM)
 - (1M stem solution: 125 µL per culture)
 - ◇ Let the cells grow further for 30min at 37°C
 - ◇ Protein expression over night at 19°C
- ◆ Cell harvest and decomposition

- ◇ Pellet the whole culture in large centrifugal jars
 - 4,500rpm, 10-20min, 4°C
- ◇ Discard the supernatant
- ◇ Resuspend the pelleted cells in 20-30mL lysis solution
 - Retransfer them into falcon tubes
- ◇ Decompose the cells using the French Press
 - 2 times for each sample
 - Pressure of approximately 1200
 - Work with ice whenever possible
- ◇ Centrifuge the samples to divide the desired substances from the cell debris
 - 4,500rpm, 1h, 4°C
 - Transfer the supernatant into falcon tubes
- ◆ Column preparation
 - ◇ Close the lower part of the columns with Para film
 - ◇ Fill them with 5mL of water(=1volume) and mark the fill level
 - ◇ Let the water out
 - ◇ Fill the column with chitin beads up to the marking
 - An evenly surface gets better results
 - ◇ Put a filter plate onto the beads
 - ◇ **From here, all steps should be carried out at 4°C**
 - ◇ Equilibrate the columns with 10volumes of column buffer (low-salt)
- ◆ Loading of the columns
 - ◇ Put the supernatant onto the columns
 - ◇ Make sure the flow-through is less than 0.5-1 mL/min

- E.g. with the use of clamps
- ◆ Washing of the columns
 - ◇ Wash the columns with 10 volumes of the column buffer (high salt)
- ◆ Induction of cleavage
 - ◇ Wash the columns with 3 volumes of cleavage buffer
 - Fast, so no clamps!
 - ◇ Let the columns induce at 4-23°C for at least 16 (up to 40) hours
- ◆ Elution
 - ◇ Eluate the columns with 1-2volumns column buffer (low salt)
 - ◇ Collect the eluate within falcon tubes or blue caps
 - ◇ Concentrate the eluate by using suitable protein filters
 - Centrifuge for 30-45minutes at full speed up to an end volume of about 250 µL
 - ◇ Wash the concentrated protein solution with 4mL protein wash buffer
 - Repeat this 2-3 times total
 - ◇ Optionally: wash with storage buffer to ensure a longer functionality
 - ◇ Make an Bradford Estimation
 - E.g. with Roti® Nanoquant([Roth](#)), see our other protocols
 - ◇ Store the proteins at the freezer or refrigerator
 - Depending on the usage
- ◆ Column regeneration
 - ◇ All of the columns can be reused 4-5 times
 - ◇ Resolve residual proteins with 3 volumes of NaOH (0.3M)
 - ◇ Incubate for 30minutes

- ◇ Wash the columns with further 7 volumes NaOH
- ◇ Wash with 20 volumes of H₂O
- ◇ Wash with 5 volumes of column buffer
- ◇ Store the columns at 4°C, closed with Para film

◆ Solutions needed for this protocol:

◇ Lysis Buffer

Substance	concentration	mass/volume in 1L
Triton X-100	0.1%	1mL
PMSF	20 μM	200 μL
TCEP	0.1-1mM	100-1000 μL
<p>Make the following stem solutions:</p> <ul style="list-style-type: none"> • 100mM PMSF • 1M TCEP <p>Resolve the substances within column buffer (low salt)</p>		

◇ Column Buffers

substance	concentration	mass/volume in 1L
Na ₂ HPO ₄	20mM	3,6g
	<ul style="list-style-type: none"> • 2 H₂O 	
NaCl	50mM (low salt)	2,92g
	500mM (high salt)	29,22g
	1000mM(high salt)	58,44g
	1500mM(high salt)	87,66g

EDTA	1mM	10mL of stem solution
<ul style="list-style-type: none"> • Prepare an EDTA stem solution of 0.1M, resolve in ddH₂O • Resolve all substances in ddH₂O • Adjust the pH to 8.0 (6.0-9.0) with phosphoric acid 		

◇ Cleavage Buffer

substance	concentration	mass/volume in 1L
Na ₂ HPO ₄	20mM	3.6g
NaCl	50mM (optionally up to 1000mM)	29g
EDTA	1mM	10mL
DTT	50mM	*
<ul style="list-style-type: none"> • *0.93g DTT resolved in 6mL ddH₂O+119mL column buffer (low salt) make up 125mL cleavage buffer • All substances resolved in H₂O • pH has to be adjusted to 7-9 with phosphoric acid • Store dark because of the DTT! 		

◇ Protein wash buffer

substance	concentration	mass/volume in 1L
Na ₂ HPO ₄	20mM	3.6g
NaCl	10mM	0.58g
EDTA	1mM	10mL

Originally from [NEB](#), now modified

