

Recombinase Polymerase Amplification (RPA) on Paper

Aim of the Experiment

Preparation of RPA mixture on paper for storage and usage in automated process on PCR Chip.

Materials

- nuclease-free H₂O (nf H₂O, Sigma Aldrich, USA)
 - 200 mM MgCl₂ (Carl Roth, Germany)
 - 280 mM MgAc (Carl Roth, Germany)
 - Liquid N₂
 - 10 μ M Primer A
 - 10 μ M Primer B
 - Template DNA
 - RNase Inhibitor, Murine (NEB, 40000 U/ml, Germany)
 - Reaction buffer (10x RNAPol Reaction Buffer, NEB, Germany)
 - T7 polymerase (T7 RNA Polymerase, NEB, Germany)
 - Ribonucleotide Solution Mix (rNTP, NEB, Germany)
 - DNaseI (RNAase-free) (NEB, Germany)
 - DNaseI buffer (NEB, Germany)
 - RPA Kit (TwistDx, UK)
 - PCR cleanup kit (NEB, Germany)
 - Polydimethylpolysiloxane (PDMS, Dow Corning, USA)
 - Heat block at 37 °C / portable PCR Device at 37 °C
 - Metal plate
-

- Dessicator
- Nuclease-free filter paper (glass microfiber paper 934-AH RTU, Whatman, GE healthcare, Germany)
- Sonicator

Procedure

1. Put metal plate into -80 °C freezer for pre-cooling. Take care that it only has little contact with the surface.
2. Mix 15 μ l of TwistDx RPA Rehydration Buffer with 2.5 μ l of Primer A and 2.5 μ l Primer B.
3. Resuspend RPA Pellet.
4. Apply resuspended solution to paperstrip that has been prepared by wax-printing to show a small circle (ca. 3 mm in diameter).
5. Carefully put the paperstrip into liquid N₂. Take care that you do not destroy the droplet that has formed on the paperstrip. Leave for several minutes.
6. Put the metal plate from the -80 °C freezer into the dessicator. Quickly apply vacuum in order to get rid of the condensed water that has crystallized on it.
7. Apply the paperstrip from the liquid N₂ to the dessicator. Leave it in the dessicator for 30 minutes.
8. Store the paperstrip in a petridish or seal in PDMS chip via sonication. Wrap with parafilm to prevent humidity to flow in.

Reaction

1. Mix 2 μ l of your sample with 29.5 μ l Rehydration Buffer, 14.5 μ l nuclease-free water and 2.5 μ l MgAc either on the PDMS chip or by pipetting. Apply to paperstrip.
 2. Incubate for 40 min at 37 °C.
 3. (Optional) For testing, use PCR clean-up kit and load on gel.
-

4. ADD GEL HERE
