

CRISPR Toxicity Assay

Overnight incubation – n=3

Inoculate a 5 mL overnight culture at 30°C of *A. tumefaciens* GV3101 with pCambiaCas.

Preparation

1. Label 30 tubes according to the following:

N1, 1 hr no plasmid	N1, 2 hr no plasmid	N1, 4 hr no plasmid	N1, 8 hr no plasmid	N1, 24 hr no plasmid
N2, 1 hr no plasmid	N2, 2 hr no plasmid	N2, 4 hr no plasmid	N2, 8 hr no plasmid	N2, 24 hr no plasmid
N3, 1 hr no plasmid	N3, 2 hr no plasmid	N3, 4 hr no plasmid	N3, 8 hr no plasmid	N3, 24 hr no plasmid
N1, 1 hr 0.5x	N1, 2 hr 0.5x	N1, 4 hr 0.5x	N1, 8 hr 0.5x	N1, 24 hr 0.5x
N2, 1 hr 0.5x	N2, 2 hr 0.5x	N2, 4 hr 0.5x	N2, 8 hr 0.5x	N2, 24 hr 0.5x
N3, 1 hr 0.5x	N3, 2 hr 0.5x	N3, 4 hr 0.5x	N3, 8 hr 0.5x	N3, 24 hr 0.5x
N1, 1 hr 1.5x	N1, 2 hr 1.5x	N1, 4 hr 1.5x	N1, 8 hr 1.5x	N1, 24 hr 1.5x
N2, 1 hr 1.5x	N2, 2 hr 1.5x	N2, 4 hr 1.5x	N2, 8 hr 1.5x	N2, 24 hr 1.5x
N3, 1 hr 1.5x	N3, 2 hr 1.5x	N3, 4 hr 1.5x	N3, 8 hr 1.5x	N3, 24 hr 1.5x

Day of

1. Dilute overnight culture 1/10 and spec it. Calculate dilutions required to achieve 0.100 abs.
2. Dilute culture according to calculation from above to achieve a total of 2 separate 35 mL cultures (minimum) that should be OD 0.100.
 - a. Culture 1: Agro grown in LB + Kan + Rif
 - b. Culture 2: Agro grown in LB + Kan + Rif + 1x Anhydrous tetracycline.
3. Split each of these cultures in to 15 sterile culture tubes, and 1 mL from the remainder of each into a cuvette.
4. Take spec reading of the remainder (this will be baseline).
5. Incubate tubes on shaker at 30°C, 220 rpm.
6. At t=2 hr, 4 hr, 8 hr, 24hr, remove designated tube from the 30°C room and spec them. Record the value.

<u>Overnight incubation</u>				Do we have?
1 colony	A. tumefaciens GV3101 w/ pCambiaCas	→	4°C room	Y
5 mL	LB media	→	Bench	Y
15 uL each	Rifampicin + Kanamycin	→	-20°C under bench	Y
1	Test tubes	→	Bench	Y

Prep				
30	Culture tubes	→	4°C Fridge	Y
Day of				
5 mL	overnight culture	→	30°C Incubator	n/a
70 mL	LB media	→	Bench	Y
65 uL	Rifampicin + Kanamycin	→	-20°C under bench	Y
35 uL	Anhydrous tetracycline stock	→	-20°C under bench	N
7, assuming we wash	Cuvettes for spec	→	Beside spec	Y
7	7 Clean cuvettes			
1 mL	LB media (labeled for CRISPR)			
	P1000			
	Sterile P1000 tips			
	Non-sterile P1000 tips			
	Pen and data sheet			
	Test tube rack			

- 12:45ish 1. Gather all the items you will need:
 7 Clean cuvettes
 LB media (labeled for CRISPR)
 P1000
 Sterile P1000 tips
 Non-sterile P1000 tips
 Pen and data sheet
 Test tube rack
- 12:50ish 2. Put pen, data sheet, 6 cuvettes and non-sterile P1000 tips by the spec machine.
- 12:55ish 3. Work near a flame and add 1 mL of LB media to a cuvette, and bring it to the spec.
 ➔ You should now have 7 cuvettes at the spec!
- 1:05ish 4. Bring test tube rack to the 30°C room and collect the 6 culture tubes, bring to spec.
- 1:05ish 5. Don't worry about being sterile – add 1 mL of each culture tube to each cuvette
 ➔ make sure you know which is which!
- 1:10pm 6. Blank the spec with the plain LB you put in a cuvette earlier. The spec should say OD600 and you shouldn't have to change any settings
 7. Read and record the OD600 value for each culture and record it on the data sheet.
 8. Bring everything back to the bench - that's it!