

Media and Agar

Aim

Prepare LB and SOC media and agar, 5x KCM, and TAE buffer.

Procedure

LB Medium

1. To make 1L LB medium, mix the following:

Table 1: LB medium mix

Component	Amount [g]
Yeast extract	5
NaCl	10
Tryptone	10
dH ₂ O	950 ml

2. Adjust the PH to 7.0 with 1 M NaOH, then autoclave.

3. Autoclave your prepared medium.

LB Agar

1. Before autoclaving, add 15 g/L agar to the LB media.
2. Cool to approximately 55 °C and add antibiotics (50 µg/µl – for 1 plate).
3. Pour into the plate (10 ml/plate), stock plates and let harden.
4. Invert and store at 4 °C.

Lab protocol

Updated: October 28th 2017
iGEM Stockholm

SOC Medium

1. To make 1L SOC medium, mix the following:

Table 2: SOC medium mix

Component	Amount [g]
Yeast extract	5
NaCl	0.585
Tryptone	20
dH ₂ O	950 ml
KCl	0.186
MgCl ₂	0.95
MgSO ₄	1.2

2. Adjust to pH 7.5 prior to use. This requires approximately 25 ml of 1 M NaOH per liter.

3. Autoclave. After cooling medium to less than 50 °C, add 20 ml filter sterilized 20 % glucose solution (4 g glucose into 20 ml).

5x KCM

1. To make 5x KCM, mix the following:

Table 3: 5x KCM Mixes

Component	Concentration [M]	Example amount [g]
KCl	0.5	37.25
CaCl ₂	0.15	16.65
MgCl ₂	0.25	50.5 (MgCl ₂ dH ₂ O)
dH ₂ O	Fill up to total volume	Fill up to 1L

2. Divide into 1 ml aliquots

1xTAE Buffer

Lab protocol

Updated: October 28th 2017
iGEM Stockholm

1. To make 50xTAE Buffer, mix the following:

Table 4: 1x TAE Mixes

Component	Concentration [M]	Example amount [g]
Tris base	-	242 g (in water)
Acetic acid	-	57.1 mL
EDTA	500 mM (pH 8.0)	100 mL
dH ₂ O	Fill up to total volume	Fill up to 1L

2. If a 1xTAE Buffer is wanted, the solution can be diluted 50:1 with dH₂O. This 1X solution will then contain 40mM Tris, 20mM acetic acid, and 1mM EDTA.

1 % Agarose mixture with GelRed dye for gel electrophoresis

1. To make a 1% Agarose mixture, mix the following:

Table 5: 1% Agarose mixture

Component	Example amount
Agarose	1g
1 x TAE Buffer	100 mL

2. After mixing the agarose and buffer in a flask according to 1, put it in the microwave and heat it for as long as it takes to completely dissolve the agarose. Be careful not to burn yourself, the liquid will become very hot. Swirl the mixture to dissolve the agarose more efficiently.

3. Add 10 μ l (1 % concentration) GelRed dye to the flask and mix it carefully.

4. Pour the agarose mixture into your casting tray and let it harden before loading your sample and running the gel electrophoresis.

Lab protocol

Updated: October 28th 2017
iGEM Stockholm

Sources

Recipes are modified from Cold Spring Harbor Protocols, 2016.

Lab protocol

Updated: October 28th 2017
iGEM Stockholm