

## Wizard® SV Gel and PCR Clean-Up System from Promega

- ◆ Gel slice and PCR product preparation
  - ◇ Dissolving the gel slice
    1. Following electrophoresis, excise DNA band from gel and place gel slice in a 1.5 mL microcentrifuge tube.
    2. Add 10  $\mu$ L **Membrane Binding Solution** per 10 mg of gel slice. Vortex and incubate at 50 – 65 °C until gel slice is completely dissolved.
  - ◇ Processing PCR amplifications
    1. Add an equal volume of **Membrane Binding Solution** to the PCR amplification.
- ◆ Binding of DNA
  - ◇ Insert **SV Minicolumn** into Collection Tube.
  - ◇ Transfer dissolved gel mixture or prepared PCR product to the Minicolumn assembly. Incubate at room temperature for 1 minute.
  - ◇ Centrifuge at 16,000 x g for 1 minute. Discard flowthrough and reinsert Minicolumn into Collection Tube.
- ◆ Washing
  - ◇ Add 700  $\mu$ L **Membrane Wash Solution** (ethanol added). Centrifuge at 16,000 x g for 1 minute. Discard flowthrough and reinsert Minicolumn into Collection Tube.
  - ◇ Add 500  $\mu$ L **Membrane Wash Solution** (ethanol added). Centrifuge at 16,000 x g for 5 minute.
  - ◇ Empty the Collection Tube and recentrifuge the column assembly for 1 minute with the microcentrifuge lid open to allow evaporation of any residual ethanol.

- ◆ Elution

- ◇ Carefully transfer Minicolumn to a clean 1.5 mL microcentrifuge tube.
- ◇ Add 50  $\mu$ L of **ddH<sub>2</sub>O** to the Minicolumn. Incubate at room temperature for 1 minute. Centrifuge at 16,000 x g for 1 minute.
- ◇ Discard Minicolumn and store at 4 °C or -20 °C.

From: [Promega](#)