

Genome editing via the CRISP-Cas9 system

Genome Editing

- ◆ Generating *E. coli* BL21 DE3 harboring pCas
 - ◇ Do a heatshock transformation of pCas with *E. coli* BL21 DE3
 - ◇ Grow the culture overnight in LB media and kanamycin at 30 °C.
 - ◇ Generate elektrocompetent *E. coli* BL21 DE3 cells harboring pCas.
 - Add Arabinose (10 mM) to the culture for lambda-Red induction.
- ◆ Transformation via electroporation of pTarget
 - ◇ Thaw 50 µL electrocompetent *E. coli* cells on ice.
 - ◇ Add 100 ng pTarget and 400 ng of the donor DNA (PCR fragment).
 - ◇ Store cells on ice for 1 minute.
 - ◇ Electroporate at $U = 2.5 \text{ kV}$, $C = 25 \text{ µF}$, $R = 400 \text{ Ω}$
 - ◇ Recover cells at 30 °C for 1 hour.
 - ◇ Plate on selective LB agar containing kanamycin (50 mg L^{-1}) and spectinomycin (50 mg L^{-1}).
 - ◇ Incubate overnight at 30 °C.
 - ◇ Verify positive transformants by colony PCR and DNA sequencing.

Plasmid curing

- ◆ Inoculate a positive colony, harboring pCas and pTarget in 2 mL LB media containing kanamycin (50 mg L^{-1}) and IPTG (0.5 mM).
- ◆ Incubate at 30 °C for 8 to 16 hours.
- ◆ Plate on selective LB agar containing kanamycin (50 mg L^{-1})
- ◆ To cure pCas, grow the cells overnight at 37 °C nonselectively.

Reference

Jiang, Y., Chen, B., Duan, C., Sun, B., Yang, J., and Yang, S. (2015). Multigene editing in the *Escherichia coli* genome via the CRISPR-Cas9 system. *Appl. Environ. Microbiol.* 81: 2506–2514.