

**SOP Name:** gBlock preparation and HiFi assembly

**Author:** Marcia Pryce

**Source(s):** Adapted from Exeter 2015

**Time Required:** 1 hour

**Notes:** HiFi assembly is equivalent to Gibson assembly except that it utilises a high fidelity polymerase

**Materials:**

- gBlocks
- Plasmid backbone (50 ng)
- TE resuspension buffer (50 µl)
- HiFi DNA assembly master mix (10 µl)
- Nuclease-free water

**Procedure:**

**gBlock preparation:**

1. Centrifuge the tube containing dried gBlocks (1000ng) for 35 seconds to ensure all DNA is at the bottom of the tube.
2. Add 50µl resuspension buffer solution for a final concentration of 20ng/µl and vortex briefly.
3. Incubate at 50C for 20 mins in a water bath.
4. Briefly vortex and centrifuge. Keep on ice, store at 20C for long term storage.

**HiFi assembly:**

1. Prepare the following assembly mix on ice:

HiFi DNA assembly Master Mix 10 µl  
Linear plasmid 50 ng  
Resuspended gBlock 2 fold increase from plasmid  
Nucleasefree water n µl  
Total volume 20µl

Use the equation below to determine the amount of gBlock required (pmols should be 2x that of the linear plasmid):  $pmols = \frac{ng\ DNA \times 1000}{bp \times 650\ Da}$

2. Incubate at 50°C for 15 minutes. Keep on ice, store at -20°C for long term storage.