

experiment log

Notebook October

zxy 昨天保存
[查看编辑历史](#) 或 [对比历史记录](#)

- 关注
- 打印
- 编辑
- 移动
- 删除

Date 10.7

Awaken of the Glycerin bacteria BL-21 containing KmAdh on the plasmid PET22b and BL-21

Recorder: Xiaoyu Zhang

Add 200 μ L of bacteria into a 5 mL LB culture and cultivate these bacteria overnight for extraction at 37 degree centigrade, 250 rpm.

Date 10.8

His-tag purification of KmADH

Recorder: Shihan Zhu

Procedure:

1. Pure out 20% ethyl inside of the tube;
2. Wash the tube with 10mL 1M imidazole 3 times;
3. Wash the tube with 10mL H₂O;
4. Balance the tube with 1xPBS buffer 3 times;
5. Load my sample;
6. Inoculate protein with Ni-NTA , 4°C, 200rpm, 1h;
7. Wash the tube with 10 mL 30mM imidazole 2 times and collect the sample;
8. Wash the tube with 10 mL 300mM imidazole and collect the sample;
9. Wash the tube with 10 mL 1M imidazole three times and collect the sample;
10. Wash the tube with 10 mL H₂O;
11. Wash the tube with 20% ethanol and seal the tube.;

PCR of KmAdh

Recorder: Xiaoyu Zhang

1.Prepare 3 PCR tubes and sequentially add :

sample	1	2	3
Sterilized ddH ₂ O	7 μ L	7 μ L	7 μ L
2 \times Primer Star	10 μ L	10 μ L	10 μ L
template	1 μ L	1 μ L	1 μ L
KmAdh-r(10 μ M)	1 μ L	1 μ L	1 μ L
KmAdh-f(10 μ M)	1 μ L	1 μ L	1 μ L
total	20 μ L	20 μ L	20 μ L

2.PCR reaction Parameters setting :

stage	temperature	time
step 1	95	10 min
step 2	98	10 s
step 3	55	15 s
step 4	72	54 s
step 5	72	10 min
step 6	4	--

30 cycles

3. Agarose gel electrophoresis mixed with X μL 6 \times DNA loading buffer each X μL sample; XX V, XX min Result:



lane1-4: Marker,WT,KmAdh,Postive control

Enzyme activity test

Recorder: Xiaoyu Zhang

Procedure: The components of enzyme active reaction system are shown below.


System components	Volume
0.1mM Tris-HCl	100 μL
3mM NADH	30 μL
40% acetaldehyde	22 μL
ddH ₂ O	873 μL
KmAdh	75 μL

The components of blank control are shown below.

System components	Volume
0.1mM Tris-HCl	100 μL
3mM NADH	30 μL
40% acetaldehyde	22 μL
ddH ₂ O	873 μL
PBS	75 μL

Scan the 340nm UV absorption value over time. Draw the OD₃₄₀-t curve. Compare the changes of OD₃₄₀ over time.

Scan the 500nm UV absorption value .



Standard curve of COBE

Date 10.15

CysDes function analysis--acid labile sulfide analysis

Recorder: Shihan Zhu


Preparation

1. Inoculate 2 mL of BL21 with CysDes into 200 mL LB media and add cysteine to final concentration 1 mM and cultivate for about 3 hours, 37°C, 250rpm.
2. When OD600 is about 0.4-0.6, add AHL to the media to final concentration 250nM.
3. Cultivate for 3 hours.

Procedure

1. 2 mL of culture samples are centrifuged at 17,000 g for 2 min;
2. The supernatant is removed and the cell pellet is resuspended in 1mL of 0.75 M NaOH;
3. The suspension is then transferred to a EP tube and incubated at 95 degree centigrade for 15 min;
4. Vortex the suspension vigorously and take 25 μ L of the suspension to mix with 375 μ L of 0.75 M NaOH and 250 μ L of 2.6% zinc acetate dihydrate ; 125 μ L 0.1% N,N-dimethyl-p-phenylenediamine dihydrochloride in 5 M HCl (freshly prepared) are added to and the solution is vortexed until clear;
5. 50 μ L of 11.5 mM FeCl₃ in 6M HCl are added and the solution is vortexed and incubated at RT for 30 min;
6. 425 μ L of deionized water are added and the OD of the samples are recorded at 670 nm. (Solutions of 0-0.2 mM sodium sulfide in 0.75 M NaOH served as calibration standards.)

SDS-PAGE



(lane left to right: before induction, after induction , raw enzyme , wt, wt)



10-20 21:59 俞文斐 创建了文档



10-22 23:19 zxy 编辑了文档 (查看更多动态)



点击发表评论

不想打字? 试试 召开视频会议

创建讨论后, 这些人会收到通知:



有新动态通知我

编辑通知成员