

Gibson Assembly

- ◆ Modified from [Gibson et al. \(2009\)](#)
- ◆ This assembly method is an isothermal, single-reaction method for assembling multiple overlapping DNA molecules. By coordinating the activity of a 5′ exonuclease, a DNA polymerase and a DNA ligase two adjacent DNA fragments with complementary terminal sequence overlaps can be joined into a covalently sealed molecule, without the use of any restriction endonuclease.
- ◆ Preparation of DNA molecules for *in vitro* recombination.
- ◆ Generate the complementary sequence overlaps by PCR using the Phusion DNA polymerase. If necessary add 5 M Betaine to the reaction mix by reducing the amount of H₂O to decrease the number of false PCR products.
- ◆ Identify the PCR products of interest by gel electrophoresis with known DNA standards.
- ◆ Extract the PCR products from the gel by cutting out the DNA fragments and clean them up by using a commercial gel clean-up kit.
- ◆ *in vitro* recombination:
- ◆ Assembly mixture:
 - ◇ 320 μL 5x isothermal reaction buffer
 - ◇ 0.64 μL of 10 U mL⁻¹ T5 exonuclease (for DNA molecules overlapping by greater than 150 bp add 3.2 μL of 10 U mL⁻¹ T5 exonuclease)
 - ◇ 20 μL of 2 U mL⁻¹ Phusion DNA polymerase
 - ◇ 160 μL of 40 U mL⁻¹ taq DNA ligase
 - ◇ add ddH₂O water up to a final volume of 1.2 mL
- ◆ Aliquot 15 μL of the reagent-enzyme mix and store it at –20 °C.
- ◆ Thaw 15 μL assembly mixture aliquot and keep it on ice until use.

- ◆ Add 5 μ l of the purified DNA molecules in equimolar amounts (between 10 and 100 ng of each DNA fragment).
- ◆ Incubate the resulting mixture at 50 °C for 15 to 60 min, with 60 min being optimal.
- ◆ Transformation (via heat shock or via electroporation) without cleaning up the assembly product.

From: [iGEM Bielefeld-CeBiTec](#)