

8-24-17 Site directed mutagenesis of Mamba in BBIC3 -> Change codon to remove 2nd EcoRI restriction site.

Purpose: Use site directed mutagenesis to change a base in the middle of EcoRI site upstream of Mamba.

Preparation

✓ Dilute mutagen primers to 125 ng/μL

✓ Dilute DNA samples to 50 ng/μL.

◦ Ligation A #4

1.2 μL sample + 8.8 μL dH₂O

◦ Ligation A #6

1.2 μL sample + 8.8 μL dH₂O

◦ Ligation A #7

1.3 μL _{sample} + 8.7 μL dH₂O

◦ Ligation A #9

2.4 μL + 7.6 μL dH₂O

8-24-17

1 Control Rxn and 4 experimental rxns.

Control Rxn volumes

5.0 NL 10X Rxn Buffer
 2.0 NL PwhiteScript DNA
 1.25 NL Control primer 1
 1.25 NL Control primer 2
 1.0 NL dNTP
 38.5 NL dH₂O
 1.0 NL Pfu Ultra HF DNA polymerase

Experimental Rxn volumes

master mix

20.0 NL 10X Rxn Buffer
 4.0 NL Antisense primer
 4.0 NL forward sense primer
 4.0 NL dNTP mix
 164.0 NL dH₂O.


Samples

4, 6, 7, 9
 from ligase A A

41.0 NL of master mix was added
 to each sample tube.

1.0 NL of each DNA experimental sample
 was added to corresponding numbered tubes.

1.0 NL of Pfu Ultra HF DNA polymerase
 was added to each tube.



8-24-17 Mamba mutagenesis.

Thermal Cycler settings

	1.	95°C	30 sec
12X Cycles	2.	95°C	30 sec
	3.	55°C	1.0 min
	4.	68°C	3.0 min

1.0 μ l of DpnI was added to the control and experimental rxns.

- centrifuged for 1.0 min.

- Incubated @ 37°C for 60 min.

DNA samples were stored in freezer overnight.

- PCR Clean up of mutant samples

- used Qiagen Qiaquick PCR clean up kit
40.0 μ l of each sample.

- kept 10.0 μ l of each sample unclean.